

Characterization of an IgG-Cleaving Protease from *Streptococcus equi* with Improved Activity Against Mouse IgGs

Chris Hosfield¹, Philip Compton², Luca Fornelli², Paul Thomas², Neil L. Kelleher², Michael Rosenblatt¹ & Marjeta Urh¹

¹ Promega Corporation, 2800 Woods Hollow Rd, Madison, WI 53711; ² Proteomics Center of Excellence, Northwestern University, Evanston, IL



1. Introduction

IdeS Protease is an immunoglobulin degrading enzyme isolated from *Streptococcus pyogenes* and it has become a valuable tool for characterization of therapeutic antibodies, Fc-fusion proteins and antibody-drug conjugates. The reasons for its value are:

- Rapid digestion (~30 minutes) of human IgGs exclusively at a single site below the hinge region yielding F(ab')₂ and Fc fragments.
- Reduction of the digestion products produces three fragments of ~25kDa that are readily analyzed by LC-MS and allow for:
 - Improved chromatographic separation of variants
 - Improved detection of modifications using high resolution MS.

Nevertheless, Ides protease shows poor activity against mouse IgGs.

IdeZ Protease was originally identified in *S. equi* subsp. *zooepidemicus*¹. Here we have expressed and purified a modified recombinant IdeZ, and show that it has significantly improved activity against mouse IgG2a and IgG3 subclasses when compared to IdeS. We also demonstrate the use of IdeZ in LC-MS workflows for human and mouse IgG characterization.

¹Lannergard & Guss, *FEMS Microbiol Lett* (2006).

Table 1. Core and lower hinge sequence of human and mouse IgG subclasses. IdeS activity against each subclass is indicated. Cleavage of human IgGs by IdeS has been previously shown to occur between two glycine residues highlighted in red.

	Subclass	Hinge/CH2 Sequence	IdeS Activity	IdeZ Activity
Human	IgG1	CPPCPAPELLGGPSVF		
	IgG2	CPPCPAPP_VAGPSVF		
	IgG3	CPRCPAPELLGGPSVF		
	IgG4	AHHAQAPEFLGGPSVF		
Mouse	IgG1	PCICTVPEV____SSVF	-	-
	IgG2a ^a	CPPCAAPNLLGGPSVF	+	++++
	IgG2b	CHKCPAPNLEGGPSVF	-	-
	IgG3	GSSCPAGNILGGPSVF	+	++++

2. Comparison of IdeS and IdeZ Activity

IdeZ & IdeS have similar performance against human and chimeric IgGs and Fc fusion proteins.

A. Therapeutic IgG Panel




Figure 1. Digestion of panel of therapeutic IgGs and Fc-fusion proteins. 50µg of therapeutic IgG or fusion protein were digested with 50 units of IdeS or IdeZ (Promega) for 30 min at 37°C in a final volume of 25µl in pH 6.6 buffer. Undigested controls and samples were analyzed by SDS-PAGE under reducing conditions.

> 95% Digestion of mouse IgG2a in 2 hours with IdeZ.

B. Mouse IgG2a Digestion






Figure 2. Digestion of recombinant mouse IgG2a with IdeS and IdeZ. 50µg of recombinant mIgG2a anti-hCD20 (Invivogen) were digested with 50 units of IdeS or IdeZ (Promega) for various time points at 37°C in a final volume of 25µl in pH 6.6 buffer. Samples were analyzed by SDS-PAGE under reducing conditions. IgG2a is completely digested in 2 hours with IdeZ whereas is only partially digested with IdeS even after O/N digestion.

3. IdeS and IdeZ Cleave Rituximab at the Same Site

IdeS and IdeZ rapidly digest IgGs below the hinge resulting in three ~25 kDa fragments after reduction.



A. TIC Chromatogram of IdeS-digested Rituximab



B. Full MS Spectra

Digestion with IdeS

