

A Complete Forensic Toxicology Screening Procedure for Drugs and Toxic Compounds in Urine and Plasma Using LC-MS/MS

Marta Kozak, Taha Rezai, Thermo Fisher Scientific, San Jose, CA

Key Words

- ToxSpec Analyzer
- ToxID Software
- LXQ Linear Ion Trap
- Clinical Toxicology
- General Unknown Screening

Introduction

Forensic toxicology laboratories commonly use automated immunoassays, gas chromatography-mass spectrometry (GC-MS) and high pressure liquid chromatography-diode array detector (HPLC-DAD) techniques to perform toxicology screening analyses. None of these techniques are able to identify all the drugs and toxic compounds that are potentially present in a sample. Implementation of liquid chromatography-mass spectrometry (LC-MS) for forensic toxicology screening provides specific and sensitive analysis of drugs and toxic substances. The benefits of the LC-MS/MS screening methodology include a simple sample preparation procedure, ease of adding new compounds to the screening method and fewer limitations based on compound volatility and thermal stability. In addition, Thermo Scientific ToxID automated toxicology screening software is able to automatically generate both Summary and Long Reports, avoiding the need for manual analysis of each sample chromatogram. This application note describes the use of the Thermo Scientific LXQ ion trap mass spectrometer equipped with an ESI source and HPLC for identification of unknown compounds in urine and plasma for clinical research and forensics.

Goal

To develop a complete LC-MS/MS forensic screening methodology which includes a sample preparation method, LC-MS method, spectra library, and data processing and reporting software.

Experimental Conditions

An MS/MS spectral library of 275 drugs and toxic compounds was created. Sample preparation of spiked human urine or human plasma was carried out using a solid-phase extraction (SPE) cartridge for basic, neutral and acidic compounds. A 13-minute LC method implementing a Perfluorophenyl (PFP) column was developed. Samples were analyzed using electrospray ionization (ESI) on an ion trap mass spectrometer in polarity switching scan dependent MS/MS experiments (see Figure 1), with retention time windows specified for each listed parent mass. The method allows acquisition of MS² spectra for co-eluting compounds and analysis of positively and negatively ionized compounds with a single run. Figure 2 shows the overall application workflow.




Figure 1: MS scan events




Figure 2: Step-by-step application workflow

Sample Preparation

Samples (1 mL of urine or 0.5 mL of plasma) were spiked with 0.1 mL of an internal standard solution at a concentration of 1 µg/mL (Chlorpromazine-D3, Haloperidol-D4 and Prazepam-D5) and diluted with 2 mL of 0.1 M phosphate buffer pH 6.0. The resulting mix was extracted with an SPE (Thermo Scientific Hypersep Verify-CX 200 mg mixed mode cartridges) procedure prior to injection onto LC-MS.

Chromatography

HPLC separation was performed with a Thermo Scientific Accela pump using a Thermo Scientific Hypersil GOLD PFP column (50 x 2.1 mm; 5 µm particles). Flow rate was set to 200 µL/min. The gradient is summarized in Table 1 (solvent A = water/0.1% formic acid/10 mM ammonium formate, solvent B = acetonitrile/0.1% formic acid). Injection volume was 10 µL.

Table 1. Thirteen-minute LC method

Time (minutes)	%A	%B
0	95	5
0.5	95	5
5.5	5	95
8.5	5	95
8.6	5	95
13	95	5

MS Conditions

Instrument:	LXQ ion trap mass spectrometer
Ionization:	ESI, Thermo Scientific Ion Max source
Capillary temperature:	275 °C
Spray voltage:	5.0 kV
Sheath gas:	30
Aux gas:	8
Data acquisition mode:	Polarity switching scan dependent experiment
Microscans:	1
WideBand Activation™:	On
Stepped Normalized	
Collision Energy:	35% ± 10%

Method Validation and Results:

The method was prequalified by processing and analyzing urine samples spiked with 10 randomly selected compounds in concentrations of 10 ng/mL, 100 ng/mL and 1000 ng/mL. Table 2 lists the concentration at which each analyte in the toxicology screen for urine samples is identified. The presence of an analyte at 10, 100 or 1000 ng/mL implies that the limit of detection is likely below that value. Of the 275 compounds analyzed, 70% were detected at 10 ng/mL, 20% at 100 ng/mL, 8% at 1000 ng/mL and 2% were detected at a concentration above 1000 ng/mL.

Table 2. Results for spiked urine samples in toxicology screen by LC-MS/MS

LXQ – 13 min method	Concentration Tested (ng/mL)		
Compound	10	100	1000
<i>All barbiturates require an APCI source for detection. P=Drug present. N=Drug not present.</i>			
11-Hydroxy-delta-9-THC	N	N	>1000
11-nor-9-carboxy-Delta-9-THC	N	N	P
2-Bromo-Alpha-Ergocryptine	P	P	P
2-Hydroxyethylflunitrazepam	N	P	P
3-Hydroxystanazolol	N	N	>1000
4-Hydroxynordiazepam	N	P	P
6-Acetylcodeine	P	P	P
6-Acetylmorphine (6-MAM)	P	P	P
7-Amino-Clonazepam	P	P	P
7-Amino-Flunitrozepam	P	P	P
Acebutolol	P	P	P
a-Hydroxy-Alprazolam	P	P	P
a-Hydroxy-Triazolam	P	P	P
Albuterol	P	P	P
alpha-Hydroxymidazolam	N	P	P
Alprazolam	P	P	P
Alprenolol	P	P	P
Aminorex	N	P	P
Amiodarone	P	P	P
Amitriptyline	P	P	P
Amlodipine	N	N	P
Amobarbital	P	P	P
Amoxapine	P	P	P
Amphetamine	P	P	P
Anhydroecgonine MethylEster	N	P	P
Antipyrine	N	N	>1000
Apomorphine	N	N	>1000

LXQ – 13 min method Compound	Concentration Tested (ng/mL)		
	10	100	1000
Astemizole	N	P	P
Atenolol	P	P	P
Atropine	N	P	P
BDB	N	P	P
Benzocaine	N	N	P
Benzoyllecgonine	N	P	P
Betaxolol	P	P	P
Bisacodyl	P	P	P
Bisoprolol	P	P	P
Bromazepam	P	P	P
Brompheniramine	P	P	P
Bupivocaine	P	P	P
Buprenorphine	P	P	P
Bupropion	P	P	P
Buspirone	P	P	P
Butalbital	N	P	P
Butorphanol	P	P	P
Cannabidiol	N	N	>1000
Cannabinol	N	N	>1000
Captopril	N	N	P
Carbamazepine	P	P	P
Carbinoxamine	N	P	P
Carisoprodol	N	N	P
Cathinone	N	N	P
Chlordiazepoxide	P	P	P
Chlorothiazide	N	P	P
Chlorpheniramine	P	P	P
Chlorpromazine	P	P	P
Chlorpromazine-D3	N	P	P
Chlorprothixene	N	N	>1000
Cinnarizine	P	P	P
cis-4-Methylaminorex	N	P	P
Cisapride	N	P	P
Citalopram	P	P	P
Clenbuterol	P	P	P
Clenbuterol	N	P	P
Clobazam	N	P	P
Clomipramine	P	P	P
Clonazepam	P	P	P
Clonidine	P	P	P
Clopidogrel	P	P	P
Clozapine	P	P	P
Cocaethylene	P	P	P
Cocaine	P	P	P
Codeine	P	P	P
Cyclobenzaprine	P	P	P
Delta9-THC	N	P	P
Desalkylflurazepam	N	P	P
Desipramine	N	P	P
Desmethyldoxepin	P	P	P
Dextromethorphan	P	P	P
Diazepam	P	P	P
Diffunisal	P	P	P
Digoxin	N	N	P
Dihydrocodeine	P	P	P
Dihydroergotamine	P	P	P
Diltiazem	P	P	P
Diphenhydramine	P	P	P
Dipyridamole	N	N	P
Disopyramide	P	P	P
Dothiepin	N	P	P
Doxepin	P	P	P
Doxylamine	P	P	P
Ecgonine-Methyl-Ester	N	N	P
EDDP	P	P	P
EMDP	P	P	P
Enalapril	P	P	P
Ephedrine	N	P	P

LXQ – 13 min method	Concentration Tested (ng/mL)			
	Compound	10	100	1000
	Ergotamine	P	P	P
	Estazolam	N	P	P
	Felcainide	P	P	P
	Fendiline	P	P	P
	Fenfluramine	P	P	P
	Fentanyl	P	P	P
	Fexofenadine	P	P	P
	Flumethasone	N	N	P
	Flunitrazepam	P	P	P
	Flunixin	N	P	P
	Fluoxetine	P	P	P
	Fluoxymesterone	N	P	P
	Fluphenazine	P	P	P
	Flurazepam	P	P	P
	Fluvoxamine	P	P	P
	Furosemide	N	P	P
	Gabapentin	N	N	P
	Gliclazide	N	N	P
	Glimepiride	N	P	P
	Glipizide	P	P	P
	Glyburide	P	P	P
	Haloperidol	P	P	P
	Haloperidol-D4	N	P	P
	Heroin	P	P	P
	HMMA	N	N	>1000
	Hydrochlorothiazide	N	N	P
	Hydrocodone	P	P	P
	Hydromorphone	P	P	P
	Hydroxyzine	N	P	P
	Imipramine	P	P	P
	Indomethacin	N	N	>1000
	Isradipine	P	P	P
	Ketamine	P	P	P
	Ketoconazole	P	P	P
	Ketoprofen	N	N	>1000
	Ketorolac	N	N	>1000
	Labetolol	N	P	P
	Lamotrigine	P	P	P
	LAMPA	P	P	P
	Lidocaine	P	P	P
	Lometazepam	N	P	P
	Loratadine	P	P	P
	Lorazepam	P	P	P
	LSD	P	P	P
	Maprotiline	P	P	P
	MBDB	N	P	P
	MDA	P	P	P
	MDEA	N	P	P
	MDMA	P	P	P
	Melatonin	N	N	>1000
	Meperidine	P	P	P
	Mepivocaine	N	P	P
	Meprobamate	N	P	P
	Mescaline	P	P	P
	Mesoridazine	P	P	P
	Metoprolol	P	P	P
	Methadienone	P	P	P
	Methadone	P	P	P
	Methamphetamine	P	P	P
	Methaqualone	N	N	>1000
	Methcathinone	N	N	P
	Methenolone	P	P	P
	Methohexital	P	P	P
	Methoxyverapmil	P	P	P
	Methylphenidate	P	P	P
	Metoclopramide	P	P	P
	Metronidazole	N	P	P
	Mexiletine	N	N	>1000

LXQ – 13 min method	
Compound	
	Mianserin
	Miconazole
	Midazolam
	Mirtazapine
	Molsidomine
	Morphine
	Morphine-3-b-glucuronide
	Nalbuphine
	Nalorphine
	Naloxone
	Naltrexone
	NAPA
	N-DemethylTrimipramine
	N-Desmethyl-cis-tramadol
	N-Desmethylflunitrazepam
	N-Desmethylselegiline
	N-DesmthylClomipramine
	N-Ethylamphetamine
	Nicardipine
	Nicotine
	Nitrazepam
	Nitrendipine
	Nizatidine
	Norbenzoylcegonine
	Norbuprenorphine
	Norclomipramine
	Norcocaethylene
	Norcocaine
	Norcodeine
	Nordiazepam
	Nordoxepin
	Norethandrolone
	Norfentanyl
	Norfluoxetine
	Norketamine
	NOR-LSD
	Normeperidine
	Normorphine
	Noroxycodone
	Noroxymorphone
	Norpropoxyphene
	Nortriptyline
	Noscapine
	OH-LSD
	Ondansetron
	Opipramol
	Oxazepam
	Oxcarbazepine
	Oxycodone
	Oxymorphone
	Papaverine
	Paraxanthine
	Paroxetine
	PCP
	Pentazocine
	Pentobarbital
	Perphenazine
	Pheniramine
	Phenobarbital
	Phenolphthalein
	Phentermine
	Phenylbutazone
	Phenyltoloxamine
	Physostigmine
	Pindolol
	Piroxicam
	PMA
	PMMA

Concentration Tested (ng/mL)		
10	100	1000
P	P	P
P	P	P
P	P	P
P	P	P
N	N	>1000
N	P	P
N	N	>1000
P	P	P
P	P	P
P	P	P
P	P	P
P	P	P
N	N	P
N	P	P
N	P	P
N	P	P
N	P	P
P	P	P
P	P	P
N	N	>1000
P	P	P
N	N	P
N	N	>1000
N	N	>1000
P	P	P
P	P	P
P	P	P
N	P	P
P	P	P
P	P	P
N	P	P
N	P	P
P	P	P
N	P	P
P	P	P
N	P	P
P	P	P
P	P	P
N	N	>1000
P	P	P
P	P	P
P	P	P
P	P	P
P	P	P
N	P	P
P	P	P
P	P	P
N	N	>1000
N	P	P
P	P	P
P	P	P
P	P	P
P	P	P
N	P	P
P	P	P
P	P	P
N	N	P
N	N	P
N	N	P
P	P	P
N	N	P
P	P	P
N	N	P
N	N	P
N	N	P
N	P	P

LXQ – 13 min method Compound	Concentration Tested (ng/mL)		
	10	100	1000
Prazepam-D5	N	P	P
Prazosin	P	P	P
Prilocaine	N	N	P
Procainamide	N	P	P
Promazine	P	P	P
Promethazine	N	P	P
Prometryn	N	P	P
Propafenone	P	P	P
Propoxyphene	P	P	P
Propranolol	P	P	P
Protriptyline	P	P	P
Psilocin	N	P	P
Pyrilamine	P	P	P
Quetiapine	P	P	P
Quinidine	P	P	P
Quinine	N	P	P
Ranitidine	N	N	P
Risperidone	P	P	P
Scopolamine	P	P	P
Secobarbital	P	P	P
Selegiline	N	P	P
Sertraline	P	P	P
Sotalol	N	P	P
Spironolactone	N	P	P
Stanozolol	N	P	P
Telmisartan	P	P	P
Temazepam	P	P	P
Terfenadine	P	P	P
Tetracine	P	P	P
Thiamylal	N	P	P
Thiopental	P	P	P
Thioridazine	P	P	P
Thiothixene	P	P	P
Timolol	P	P	P
Topiramate	P	P	P
Trazodone	P	P	P
Triazolam	P	P	P
Trimethoprim	P	P	P
Trimipramine	P	P	P
Venlafaxine	P	P	P
Verapamil	P	P	P
Vincristine	P	P	P
Warfarin	P	P	P
Zimelidine	P	P	P
Zolpidem	P	P	P
Zopiclone	N	N	P

All barbiturates require an APCI source for detection. P=Drug present. N=Drug not present.

Table 3. Results for spiked plasma samples in toxicology screen by LC-MS/MS

LXQ – 13 min method Compound	Concentration Tested (ng/mL)		
	10	100	1000
BDB	N	P	P
Benzocaine	N	P	P
Benzoyllecgonine	P	P	P
Betaxolol	P	P	P
Bisacodyl	P	P	P
Bisoprolol	P	P	P
Bromazepam	N	P	P
Brompheniramine	N	P	P
Bufotenine	N	P	P
Bupivocaine	P	P	P
Buprenorphine	P	P	P
Bupropion	N	P	P
Buspirone	P	P	P
Butorphanol	P	P	P
Cannabidiol	N	P	P
Cannabinol	N	P	P
Captopril	N	N	>1000
Estazolam	N	P	P
Carbamazepine	P	P	P
Carbinoxamine	P	P	P
Carisoprodol	N	P	P
Cathinone	N	N	>1000
Chlordiazepoxide	N	P	P
Chloroquine	N	P	P
Chlorpheniramine	P	P	P
Chlorpromazine	N	P	P
Chlorprotixene	P	P	P
Clozapine N-Oxide	N	P	P

All barbiturates require an APCI source for detection. P=Drug present. N=Drug not present.

For selected sets of compounds the method was also prequalified by processing and analyzing spiked plasma samples. Table 3 lists the concentration at which each analyte in the toxicology screen for plasma samples is identified. In general, detection limits for urine and plasma are comparable.

In addition, the assay performance was verified by analyzing urine samples and data were compared to the results from established LC-UV and immunoassay analytical techniques. The result is shown in Table 4. The LC-MS/MS method has consistently identified more analytes present in the sample than either LC-UV or immunoassays.

Table 4. Urine sample analyzed with LC-MS/MS, LC-UV and Immunoassay methods

LC-MS	LC-UV	Immunoassay
Nortriptyline	Nortriptyline	Barbiturates
Amitriptyline	Amitriptyline	Benzodiazepines
Benzoyllecgonine	Benzoyllecgonine	Cocaine
Cocaine	Cocaine	Opiates
Norcoethylen	Cocaehtylene	THC
Norbenzoyllecgonine	-	-
Morphine	-	-
Norcocaine	-	-





Figure 3: The ToxID Summary Report

Company Name
ToxID Summary Report

C:\Documents and Settings\marta.kozak\Desktop\Desktop\Application_Notes\ToxID\2J.RAW
C:\Documents and Settings\marta.kozak\Desktop\Desktop\Application_Notes\ToxID\ToxID_config_13min.csv



Code	SI	RSI	m/z	Expected RT	Real RT	Intensity	Library Name
p	909	909	240.0	5.20	5.06	148721	Tox_Library
p	857	873	278.2	6.60	6.76	9549	Tox_Library
p	816	837	278.2	4.90	4.83	12964	Tox_Library
p	932	932	310.2	6.70	6.83	42262	Tox_Library
i	859	859	322.2	6.80	6.82	2924	Tox_Library
i	969	974	330.2	5.60	5.56	154827	Tox_Library
i	830	837	380.4	6.20	6.26	85589	Tox_Library
p	870	871	384.2	5.40	5.45	39512	Tox_Library

port is designed for a quick synopsis of the data.


Table 5. Simple workflow for adding new analytes

STEP 1: Directly infuse analyte to obtain MS ² spectra, then add spectra to the library	10 Minutes
STEP 2: Run analyte on column to obtain retention times	13 Minutes
STEP 3: Update Parent Mass Table in instrument method with parent masses and retention times	2 Minutes
STEP 4: Update ToxID with name, parent masses, the most intense product ion and retention times	2 Minutes


Company Name ToxID Long Report

Raw File Name: C:\Documents and Settings\marta.kozak\Desktop\Desktop\Application_Notes\ToxID\2J.RAW
Config File Name: C:\Xcalibur\examples\ToxID\ToxID_config_13min.csv
Sample Name:
Laboratory: ChemLab
Acquisition Start Time: 2/13/2007 1:04:54 AM


Peak Number	Compound Name	Code	SI	RSI	m/z	Expected RT	Real RT	Intensity	Library Name
3	Venlafaxine	p	816	837	278.2	4.90	4.83	12964	Tox_Library




Acquired Spectrum



Delta Spectrum



Library Spectrum



Library Structure




Figure 4: The ToxID Long Report is designed for a more thorough examination of the data.

ToxID™ Software Automates Reporting, Reduces Manual Analysis

ToxID software identifies compounds present in the sample based on MS/MS spectra and retention times. Positive hits are automatically reported via ToxID software. Reports are automatically generated, reducing the time necessary for manual analysis of each sample chromatogram. An example of a Summary Report is shown in Figure 3. A Long Report with one page per detected compound is shown in Figure 4.

Adding New Compounds to the Application

This LC-MS/MS workflow allows the user to quickly and easily add new analytes to the screening method. This feature is very important for forensic toxicology screening because new target compounds are continually being added to the target list. As shown in Table 5, new compounds can typically be added in less than 1 hour.

Conclusion

The comprehensive, turn-key forensic toxicology screening methodology described in this application note utilizes an LXQ ion trap, and includes an SPE procedure and LC method that enables the identification of 275 compounds in human urine and human plasma. Accompanying ToxID software performs automatic data analysis and reporting. This eliminates the need for manual data interpretation and increases confidence in compound identification. It is worth noting that when compared to other screening methods, the LC-MS/MS screening methodology identifies more analytes.

For Research Use Only. Not for use in diagnostic procedures.

Legal Notices

©2009 Thermo Fisher Scientific Inc. All rights reserved. All trademarks are the property of Thermo Fisher Scientific Inc. and its subsidiaries. This information is presented as an example of the capabilities of Thermo Fisher Scientific Inc. products. It is not intended to encourage use of these products in any manners that might infringe the intellectual property rights of others. Specifications, terms and pricing are subject to change. Not all products are available in all countries. Please consult your local sales representative for details.

View additional Thermo Scientific LC/MS application notes at: www.thermo.com/appnotes

In addition to these offices, Thermo Fisher Scientific maintains a network of representative organizations throughout the world.

Africa-Other

+27 11 570 1840

Australia

+61 2 8844 9500

Austria

+43 1 333 50 34 0

Belgium

+32 2 482 30 30

Canada

+1 800 530 8447

China

+86 10 8419 3588

Denmark

+45 70 23 62 60

Europe-Other

+43 1 333 50 34 0

Finland/Norway/Sweden

+46 8 556 468 00

France

+33 1 60 92 48 00

Germany

+49 6103 408 1014

India

+91 22 6742 9434

Italy

+39 02 950 591

Japan

+81 45 453 9100

Latin America

+1 608 276 5659

Middle East

+43 1 333 50 34 0

Netherlands

+31 76 579 55 55

South Africa

+27 11 570 1840

Spain

+34 914 845 965

Switzerland

+41 61 716 77 00

UK

+44 1442 233555

USA

+1 800 532 4752

www.thermo.com



Thermo Fisher Scientific,
San Jose, CA USA is ISO Certified.

AN62924_E 12/09S