

# ADP-Glo™ Kinase Profiling Systems for targeted and flexible kinase inhibitor profiling

Jacquelyn Hennek, Juliano Alves, Said Goueli, and Hicham Zegzouti

Promega Corporation, 2800 Woods Hollow Road, Madison, WI 53711





## 1. Abstract

Drug safety is of paramount importance in the pharmaceutical industry indicating that minimal side effects constitute a major requirement in drug development. Therefore, novel drug candidates need to be profiled against various liability targets, including a broad panel of kinases to provide a better understanding of off-target activities. Potentially, profiling can also identify new targets that may lead to novel therapeutic indications. A universal, robust and affordable technology is desirable to assess selectivity and potency of drug candidates against multiple classes of kinases. The luminescent ADP-Glo™ kinase assay is a universal platform that measures kinase activity by quantifying the amount of ADP produced during the enzymatic reaction. We have tested the utility of this platform with 174 optimized Kinase Enzyme Systems (KES) spanning different families of the human kinome. Here we present standardized Kinase Profiling Systems for simple kinase inhibitor profiling studies. The Kinase Profiling Systems are a set of kinases organized by families and presented in easy to use multi-well strips. Each strip contains eight enzymes each with their corresponding substrates, and standardized for optimal kinase activity for inhibitor profiling. Using the profiling strips we easily generated selectivity profiles, identifying compound promiscuity towards members of a single kinase subfamily or different subfamilies of the kinome. The ADP-Glo™ KES platforms now address the needs of basic kinase characterization, kinase screening, mode of action (MOA) studies and profiling in an affordable manner using one assay format.

## 2. ADP-Glo™ is a positive detection Assay for product formation.

### Assay concept, formats and Features




- Universal:** Any kinase-substrate combination and ATPases.
- Wide dynamic range:** High sensitivity at low % ATP to ADP conversion allows use of lower amount of enzyme.
- Broad range of [ATP]: ( $\mu$ M to mM)** allows distinction between ATP competitive and non competitive inhibitors.
- Robust Assay:** Homogenous with Z' higher than 0.7.

## 3. Promega validated kinase panel covers the human Kinome.

Kinase Enzyme System (KES) is a complete Kinase assay solution

ADP-Glo™ Kinase Assay 0-1mM ATP	
0.5ml	10mM UltraPure ATP
0.5ml	10mM ADP
5ml	ADP-Glo Reagent
10ml	Kinase Detection Buffer
1cake	Kinase Detection Substrate


Akt1 Kinase Enzyme System (Example)	
0.1ml	Akt1 Kinase (10 $\mu$ g)
1ml	AKT (SGK) substrate (1mg)
1.5ml	Kinase Assay buffer
25 $\mu$ l	100mM DTT
25 $\mu$ l	2.5mM MnCl <sub>2</sub>
500 $\mu$ l	Kinase Activator**



Broad Human Kinome coverage with >170 KES

## 4. Optimal performance ensured by reagent purity.

Genuine Kinase activities validated through specific inhibitions and intrinsic ATPase activity evaluation




- Kinases have different levels of intrinsic ATPase activity in vitro.
- The Kinases and their ATPase activities are completely inhibited by pan or selective inhibitors consistent with high purity.
- Specific inhibitors have the same potency (IC<sub>50</sub>) towards the kinase or its intrinsic ATPase activity confirming the accuracy of inhibitor profiles.

## 5. Kinase strips make profiling with ADP-Glo™ platform simple.

Important kinase targets organized in multi-well strip panels


1	2	3	4	5	6	7	8	9	10	11	12	General Panel
EGFR	ABL1	AKT1	CDK2/CyclinA2	ERK2	AKT1	CHK1	Aurora A	AMPK <sub>α1</sub>	AKT2	ASK1	FGR1	CDK2/CyclinA2
HER2	BRK	EPHA1	CDK2/CyclinA2	HER2	CDK2/CyclinA2	CHK2	Aurora B	AMPK <sub>α2</sub>	CDK2/CyclinA2	HPK1	JAK3	CDK2/CyclinA2
HER4	BTK	FAK	CDK2/CyclinA2	ITK	CDK2/CyclinA2	CDK5/p35	CDK5/p35	CDK6/CyclinD3	CDK6/CyclinD3	CDK9/CyclinK	PKC <sub>α</sub>	CDK2/CyclinA2
IGF1R	CSK	CDK5/p35	CDK5/p35	JNK1	CDK5/p35	CDK5/p35	CDK5/p35	CDK9/CyclinK	CDK9/CyclinK	CK2 <sub>α</sub>	PKC <sub>β</sub>	CDK2/CyclinA2
InR	FYN A	JAK3	CDK5/p35	p38 $\alpha$	CDK5/p35	CDK5/p35	CDK5/p35	CDK9/CyclinK	CDK9/CyclinK	CK2 <sub>β</sub>	PKC <sub>γ</sub>	CDK2/CyclinA2
KDR	LCK	PYK2	CDK5/p35	p38 $\beta$	CDK5/p35	CDK5/p35	CDK5/p35	CDK9/CyclinK	CDK9/CyclinK	CK2 <sub>α</sub>	PKC <sub>δ</sub>	CDK2/CyclinA2
PDGFR $\alpha$	LYN B	SYK	CDK5/p35	p38 $\delta$	CDK5/p35	CDK5/p35	CDK5/p35	CDK9/CyclinK	CDK9/CyclinK	CK2 <sub>β</sub>	PKC <sub>ε</sub>	CDK2/CyclinA2
PDGFR $\beta$	SRC	TKA	CDK5/p35	RSK2	CDK5/p35	CDK5/p35	CDK5/p35	CDK9/CyclinK	CDK9/CyclinK	CK2 <sub>α</sub>	PKC <sub>η</sub>	CDK2/CyclinA2

### Streamlined profiling protocol



## 6. Simple protocol for flexible and targeted inhibitor profiling.

Dose response for 2 compounds against 1 Kinase Profiling Strip



One Kinase Profiling Strip tested with 8 compounds at a single dose

Kinase	Gefitinib	Dasatinib	Tofacitinib	SB20350	Roscovitine	PF-477736	Tozertib	Enzastaurin
AKT1	108	92	105	94	102	54	107	101
EPHA1	84	1	102	79	97	63	77	100
Fak	95	84	118	106	109	86	91	115
ITK	102	101	99	94	100	33	24	99
JAK3	99	113	0	87	94	36	88	103
PYK2	133	134	122	106	134	120	80	116
SYK	104	72	84	87	94	34	77	103
TRKA	95	82	60	89	87	9	1	89

Inhibitor Concentration = 1  $\mu$ M

- Enabling flexible Kinase Profiling with the Strip Systems.
- Dose response or single dose profiling against 8 kinases at once.

## 7. Creating selectivity profiles of inhibitors with ADP-Glo™ platform.

Inhibitor profiling performed on a large scale using a single dose compound profiling protocol

Strip #	Kinase	ADP-Glo		Enzastaurin		Gefitinib		Roscovitine		Tozertib		Kinase	
		ST-1	ST-2	ST-3	ST-4	ST-5	ST-6	ST-7	ST-8	ST-9	ST-10	ST-11	ST-12
EGFR	99	104	43	4	64	102	94	103	21	96	97	78	EGFR
HER2	102	106	102	102	104	101	95	103	101	103	102	102	HER2
HER4	95	104	26	26	97	102	97	97	97	98	97	97	HER4
IGF1R	102	100	98	103	100	102	97	97	97	98	97	97	IGF1R
InR	98	100	95	97	97	97	97	97	97	97	97	97	InR
PDGFR $\alpha$	84	93	86	86	84	83	81	84	84	84	84	84	PDGFR $\alpha$
PDGFR $\beta$	85	88	85	85	87	87	87	87	85	85	85	85	PDGFR $\beta$
ALK1	98	102	27	97	97	106	107	20	97	97	97	97	ALK1
BRK	95	105	36	73	95	101	95	103	101	103	101	101	BRK
CSK	95	103	41	84	93	99	92	81	99	98	97	97	CSK
FYN													