

Frauke Henjes<sup>1</sup>, Lise Aagaard<sup>1</sup>, Darren Ellis<sup>1</sup>, Hannah Lindstroem<sup>1</sup>, Morten Oksvold<sup>2,3</sup>, Axl Neurauter<sup>1</sup>, Ketil W. Pedersen<sup>1</sup>.

<sup>1</sup>Life Technologies AS, Ullerhausséen 52, 0379, Oslo, Norway; <sup>2</sup>Department of Immunology, Institute for Cancer Research, Oslo University Hospital, Oslo, Norway;

<sup>3</sup>Centre for Cancer Biomedicine, University of Oslo, Oslo, Norway


Introduction

Invasive biopsies could often be replaced by simpler and safer non-invasive liquid biopsies using plasma or blood samples. A non-invasive liquid biopsy permits the analysis of multiple circulating biomarkers in one sample and facilitates the discovery of disease as well as treatment response monitoring in applications such as oncology. Here, we demonstrate that paramagnetic Dynabeads® provide a versatile and automation friendly tool for fast, sensitive, and high throughput isolation that can be tailored to isolate specific circulating biomarkers such as cells, exosomes, and cell-free nucleic acids.

Circulating Tumor Cells (CTCs)

Positive Isolation


For the positive isolation, Dynabeads® coupled with a monoclonal antibody targeting EpCAM (epithelial cell adhesion molecule) are used to capture CTCs from whole blood samples. For downstream applications, mRNA can be isolated using Dynabeads® Oligo (dT)<sub>25</sub> and analyzed using gene expression analysis (Figure 1E) or sequencing.





Negative Isolation

CTCs can also be enriched by applying negative isolation. The advantage is that it is independent of surface marker expression (e.g. EpCAM) on the CTCs. To deplete white blood cells, Dynabeads® CD45 are used. These beads deplete more than 95 % of the leukocytes from mononuclear cells (MNCs) (Figure 2), whole blood, and buffy coat. Addition of Dynabeads® CD15 could further increase cell depletion. Following negative isolation, downstream applications could include cell culturing and molecular analysis such as qPCR and sequencing.

**Figure 2:** Fluorescence histograms for CD45-stained MNCs with inlaid scatter plots, measured before **(A)** and after **(B)** depletion of CD45-positive cells using Dynabeads® CD45.



Exosomes




**Figure 4:** **A)** Transmission Electron Microscopy (TEM) of pre-enriched Sudhl4 exosomes labeled for CD63 <sup>2)</sup>. CD63 was detected using anti-CD63 antibodies followed by protein A gold (10 nm). **B)** TEM of exosomes isolated with Dynabeads® Human CD63 Isolation/Detection (human cell culture media)<sup>1)</sup>. Bar 500 nm.

Circulating Nucleic Acids

In order to isolate sufficient amounts of cfDNA for downstream analysis, typically 4 to 10 mL plasma or serum is necessary to be processed into volumes below 100 µL. During processing, the work volumes can increase over 3 fold to that of the input volume. Processing is therefore mostly done manually and is highly work and skill intensive, while automation presents a serious challenge. Here we provide a solution using magnetic beads, Dynabeads® MyOne™ SILANE, in a highly sensitive, high-throughput automatable large volume sample processing.

**Figure 5:** Dynabeads® MyOne™ SILANE Automation on Tecan® Freedom Evo® 150



**Figure 7:** Four mL plasma samples from 10 individuals (A-J) and pooled plasma derived from Cell-Free DNA BCT tubes (S11 and 12) processed manually and on an automated setting (Figure 5). The resulting eluates were then quantified by qPCR with fragment specific probes.

**A/B)** Recovery rates of samples laced with 10<sup>6</sup> copies of 120 bp **(A)** and 220 bp DNA **(B)** (~100%). **C)** Recovery rates of endogenous cfDNA of manual and automated sample processing in comparison to that of a non-automatable Competitor. Dynabeads® MyOne™ SILANE yielded equally efficient recovery manually and automated, albeit yielding slightly lower recovery rates than that of Competitor (87% and 83% manually and on Tecan® instrument, respectively, to Competitor). **D)** As a guideline for sample eluate purity, the total protein contents of Dynabeads® MyOne™ SILANE manual and automated extractions and Competitor were compared. The lowest mean total protein content was achieved with Dynabeads® MyOne™ SILANE manual processing. Competitor gave medium protein content but showed high sample to sample variation.

**Figure 6:** SILANE cfDNA Isolation from 4 mL Plasma or Serum

