



Human iPSC-derived cardiomyocytes: A translational model to predict drug-induced cardiac arrhythmias and long-term toxicity

R Kettenhofen¹, G Luerman¹, R Bucerius¹, I Kopljär², DJ Gallacher², A De Bondt², E Vlamincx², I Van Den Wyngaert², HR Lu², H Bohlen¹

¹ Axiogenesis AG, Cologne, Germany; ² Discovery Sciences, Janssen Research and Development, Janssen Pharmaceutica, Beerse, Belgium
Conference Contact: roman.bucerius@axiogenesis.com

BACKGROUND


Human induced pluripotent stem cell (iPSC)-derived cardiomyocytes represent a promising model for *in vitro* prediction of cardiac arrhythmias. Currently, two commercially available hiPSC-derived cardiomyocyte products - including Axiogenesis Cor.4U human cardiomyocytes - are validated in an international multi-coresite study inside the Comprehensive *in vitro* Pro-arrhythmia (CiPA) consortium. This aims to change the paradigm of safety pharmacological assessment from the assessment of drug interaction with the hERG channel to a potentially clinically more relevant assay system using human iPSC-derived cardiomyocytes in MEA-, voltage- and calcium-sensitive dye recordings.

Besides CiPA, a variety of oncological drugs including tyrosine kinase and HDAC inhibitors have been reported to induce long-term cardiac toxicity in the clinic. Recent findings clearly show that Cor.4U cardiomyocytes can be applied to a long-term impedance assay and the results from known drugs perfectly translate to the clinical observations.


The results shown here imply that iPSC-derived Cor.4U human cardiomyocytes are a translational *in vitro* cell model for the prediction of clinically relevant drug-induced cardiac arrhythmias and long-term toxicity.

RESULTS


RECENT PUBLICATIONS: ACUTE ION CHANNEL PHARMACOLOGY AND ARRHYTHMIA ASSESSMENT



FUNCTIONAL AND TRANSCRIPTIONAL CHARACTERIZATION OF HISTONE DEACETYLASE INHIBITOR-MEDIATED CARDIAC ADVERSE EFFECTS IN COR.4U CARDIOMYOCYTES



CHRONIC CARDIAC TOXICITY ASSOCIATED WITH TKIs - ASSESSMENT WITH CARDIOMYOCYTE IMPEDANCE ASSAY




CONCLUSIONS

- Human iPSC-derived Cor.4U cardiomyocytes can be employed in a number of higher throughput phenotypic screening assays (voltage dyes, microelectrode arrays, impedance) to identify acute ion channel/arrhythmia concerns, chronic cardiotox effects from a number of mechanisms (HDAC, TKI, structural concerns, etc).
- Cor.4U combined with the impedance system allow for long-term as well as label-free cardiotoxicity assays on a medium throughput level.
- HDAC inhibitors can be triaged / prioritized according to their *in vitro* cardiotoxicity profiles as the phenotypic impedance assay/beating characteristics and transcriptional / genetic changes exhibit responses consistent with clinical (cardiotox) observations.
- Overall, human iPSC-derived cardiomyocytes present an ideal tool for acute as well as chronic assessment of drug-induced cardiotoxicity while new optogenetics tools and integrated sensors will usher in a new era of high throughput cardiotoxicity assessment.

FOR MORE INFORMATION VISIT WWW.AXIOGENESIS.COM OR CONTACT INFO@AXIOGENESIS.COM


NEW TOOLS TO ACCELERATE DRUG DISCOVERY

Optical pacing of Chr2-YFP expressing Cor.4U on Nanion CardioExcyte 96



Simultaneous recording of impedance and extracellular field potential (EFP) during optical stimulation. Cells were stimulated on CardioExcyte NSP96 plates applying 2 ms blue light pulses. Optical stimulation allows for overlay of the signal giving in-depth mechanistic understanding of the contractility and electrophysiology of iPSC-derived cardiomyocytes.

Recording of calcium transients in GCaMP6f expressing Cor.4U on Hamamatsu FDSS 7000EX



Measurement of GCaMP6f fluorescence intensity during changes in intracellular calcium concentrations in Cor.4U cells. Left: Simultaneous acquisition of fluorescence intensity in all 96 wells. Nearly all wells showed low background fluorescence and high signal to noise. The duration of calcium transients varied between 550 - 600 ms with beat rates comparable to untreated cells (70 bpm). Right: Magnified fluorescence intensity profile. Images 1 - 5 demonstrate high transfection efficiency of GCaMP6f and fluorescence signal during cell contraction.