

Analysis of Urine SRMs Using Solid-Phase Micro Extraction, Dynamic Headspace, and Liquid Injection with Comprehensive Two-Dimensional Gas Chromatography (GCxGC) High Resolution Time-of-Flight Mass Spectrometry

David E. Alonso, Joseph E. Binkley, and Jonathan D. Byer | LECO Corporation, St. Joseph, MI

Introduction

Introduction

Laboratory medicine began thousands of years ago through urine analysis. Urine was a "Divine Fluid and Window to the Body". During Babylonian, Egyptian, and through Victorian times, urine was the primary diagnostic tool (Uroscopy).

Today urine is still a favored biofluid for medical diagnostic testing (Urinalysis) because large volumes are easily obtained. In addition, urine is relatively free from interfering proteins and lipids, and it tends to "hold" high concentrations of drugs and metabolites over extended periods of time. Modern, routine clinical tests include the determination of specific gravity and measurement of levels of glucose, nitrates, etc.

Objectives

- o Evaluate different methods of sample introduction
- o Implement the use of enhanced, comprehensive (GCxGC) chromatography for separation of urine compounds
- o Use modern-day, high resolution time-of-flight mass spectrometry (HRT) and powerful processing software to quickly and confidently identify compounds in urine

Technology

Sample Introduction



Data Acquisition & Processing

LECO Pegasus® GC-HRT & HRT 4D with ChromaTOF-HRT® Brand Software



Experimental

Sample Preparation

- a) Solid Phase Micro Extraction (SPME), Twister, and Dynamic Head Space (DHS) analysis:

1 mL aliquots of smoker's (NIST SRM 3673) and non-smoker's (SRM 3672) urine were placed in 10 mL vials. Salt (200 mg NaCl) was added to increase the ionic strength of the mixture. The urine was sampled without further treatment or after the addition of either acid (60 µL 12M HCl) or base (100 µL 5M NaOH).

- b) Liquid sampling:

200 µL of urine incubated with 10 mg of urease at 37°C for 15 minutes, 800 µL of methanol was added, and the mixture was vortexed and then centrifuged (12,000 g for 10 minutes). 500 µL aliquots were evaporated to dryness and derivatized using a two-step process:

- 1) Methoximation (30 µL of MEOX, 60°C, 1 hr.) and
- 2) Silylation (100 µL of ACN, 200 µL MSTFA, 60°C, 1 hr).

Instrument Parameters

- **SPME:** Extracted samples with an 85 µm polyacrylate-coated SPME fiber for 30 minutes in an agitator at 80°C (5 s on, 2 s off, 250 rpm). Desorbed for 3 minutes into a CIS4 injector at 280°C.
- **DHS:** Sampled urine for 1 hour using Tenax® desorption tubes (complete evaporation). Desorbed compounds using a TDU (40°C to 280°C at 720°C/min, 5 min hold), and CIS4 injector (-120°C to 280°C at 720°C/min, 5 min hold).
- **Twister:** Extracted samples with PDMS coated stir bars for 1 hour and transferred bars to TDU tubes. TDU Desorption (30°C to 280°C at 12°C/s, 5 min hold), and CIS4 injector (-30°C to 280°C at 12°C/s, 5 min hold).
- **Liquid:** 1 µL injection, split 20:1, CIS4 injector (35°C to 280°C at 12°C/s, 5 min hold)

Pegasus GC-HRT 4D

| Gas Chromatograph | Agilent 7890B with Gerstel MPS Autosampler |
|-------------------------------|--|
| Injection | 1µL, Split 20:1, CIS 4 (35 °C to 280 °C at 12 °C/s, hold 5 minutes) |
| Carrier Gas | He @ 1.0 mL/min, Constant Flow |
| Column 1 | Rxi-5 Sil MS, 30 m x 0.25 mm i.d. x 0.25 µm (Restek, Bellefonte, PA, USA) |
| Column 2 (GCxGC-HRT) | Rxi-17 Sil MS, 0.6 m x 0.25 mm i.d. x 0.25 µm (Restek, Bellefonte, PA, USA) |
| Temperature Program | Liq.: 50 °C (1 min) to 150 °C @ 10 °C/min, to 320 °C @ 15 °C/min (10 min) SPME,DHS,Twister: 50 °C (1 min) to 300 °C @ 10 °C/min, to 320 °C @ 30 °C/min (6 min) |
| Modulation (GCxGC-HRT) | 3s with 5 °C secondary oven temperature offset and 15 °C modulator offset |
| Mass Spectrometer | LECO Pegasus® HRT |
| Transfer Line | 300 °C |
| Ion Source Temperature | 250 °C (EI) |
| Acquisition Mode | High Resolution, R = 25,000 at m/z = 219 |
| Ionization & Mass Range (m/z) | EI35-510 |
| Acquisition Rate | 10 spectra/s; (200 spectra/s GCxGC-HRT) |

Results: SPME




Fig. 1. AIC and Table of Representative Compounds in Smoker's Urine.




Fig. 2. Peak True (Deconvoluted) and Library Mass Spectral Data for p-Cresol and Androst-2-en-17-one.




Fig. 3. Peak True and Library Mass Spectral Data for L-Borneol and α-Corocalene.

Results: DHS




Fig. 4. AIC and Table of Representative Compounds in Smoker's Urine.




Fig. 5. AIC and XIC Expansion of Smoker's Urine Demonstrating benefits of High Resolution Deconvolution. Peak True and Library Mass Spectral Data for Palmitic Acid and Dibutyl Phthalate.

Results: Twister




Fig. 6. AIC and Table of Representative Compounds in Smoker's Urine. Example of HRD showing successful deconvolution and identification of Dehydroandrosterone Acetate and TDCPP.

Results: Derivatization/Liquid Injection

Pegasus GC-HRT 4D: GCxGC-HRTOFMS

- ✓ Enhanced Chromatographic and Mass Spectral Resolution
- ✓ Group Clustering – Structured Chromatograms
- ✓ Improved Characterization of Compounds




Fig. 7. GCxGC-HRT OFMS Instrument.




Fig. 8. A) GC-HRT AIC and B) GCxGC-HRT 4D Contour Plot Smoker's Urine.




Fig. 9. A) Peak True; B) Library Mass Spectral Data for Gabapentin/Lactam in Smoker's Urine; C) Table of Representative Compounds; D) Contour Plot.




Fig. 10. Contour Plot of Non-Smoker's Urine. Table of Representative Mono-Nucleosides and Disaccharides in Non-Smoker's Urine.

Conclusion

- Enhanced two-dimensional chromatographic resolution, spectral similarity searches of large, well-established databases, and formula determinations using high resolution accurate mass ions increases confidence in metabolite identification.
- The Pegasus GC-HRT and Pegasus GC-HRT 4D facilitate fast and confident compound identification by utilizing industry leading high resolution deconvolution (HRD).
- SPME, DHS, and Twister compound sampling techniques are viable alternatives for the analysis of volatiles, semi-volatiles, and hydrophilic compounds in urine samples.
- Derivatization and liquid injections greatly increased the number of compound classes amenable to GC-HRT and GC-HRT 4D analyses.