

Improved Ligation Specificity with Chemically Modified Ligation Components

Sabrina Shore, Alexandre Lebedev, Elena Hidalgo Ashrafi, Gerald Zon, Natasha Paul, Richard Hogrefe.

TriLink BioTechnologies, Inc. San Diego, CA 92121.

Abstract

Ligases are gaining utility in molecular biology applications, such as nucleotide sequence detection, single nucleotide polymorphism (SNP) detection, protein detection and "next generation" sequencing by ligation. With the increased demand for DNA ligases in the field of biotechnology, comes increased demand for ligation fidelity. Described approaches to improved ligation fidelity include ligases from different biological sources, point mutations of key amino acid residues within the ligase, modified reaction conditions and addition of crowding reagents, such as PEG. Although most approaches to improved ligation fidelity have focused on the ligase itself, further improvements are needed and may be attainable by a different approach. Herein a strategy to improve the discrimination between matched and mismatched targets is described which employs chemical modification to the nucleic acid components of the reaction, such as the donor probe, the acceptor probe and the ATP cofactor. The results demonstrate that chemically modified components increase the stringency of DNA ligase-mediated nucleic acid detection, providing a unique approach for SNP genotyping.

Figure 1

Proposed approach for improved DNA ligase specificity




Figure 2

Limitations of T4 DNA ligase in mismatch discrimination




Figure 3

Design of chemically modified ATP Cofactors





Figure 4


Evaluation of chemically modified cofactor analogs



Ligation conditions: 1X Ligase Buffer (50 mM Tris (pH 7.5), 10 mM DTT, 10 mM MgCl₂), Ligation Probes (1 μM), ssDNA Template (1 μM), ATP analogs (1 mM), T4 DNA Ligase (0.5 Weiss units), 20 μL. Thermal cycling conditions: 95°C (3 min), 4°C (3 min), add Ligase, 22°C (20 min), 65°C (10 min).

Figure 8


Evaluation of lead probe and cofactor modifications with T4 DNA Ligase



Ligation conditions: 1X T4 Ligase buffer (50 mM Tris (pH 7.5), 10 mM DTT, 10 mM MgCl₂), Ligation probes (1 μM), ATP analogs (1 mM), ssDNA template (1 μM), T4 DNA Ligase (0.5 Weiss units), 20 μL. Thermal cycling conditions: 22°C (1 hour), 65°C (10 min).

Figure 9


Real-time ligation-PCR using chemically modified ligation components



PCR conditions: 1X PCR buffer (20 mM Tris (pH 8.4), 50 mM KCl, 2.5 mM MgCl₂), CleanAmp™ Precision Primers (0.1 μM), 0.2 mM dNTPs, 5 μL Ligation product (diluted 10³), 1.25 μL iVGN Taq DNA polymerase, SYTO®9 Green (2 μM), ROX (0.03 μM), 25 μL. Thermal cycling conditions: 95°C (10 min); [95°C (30 sec), 56°C (30 sec), 72°C (1 min)] 40X.

Figure 10

Cystic Fibrosis SNP detection by Ligase Detection Reaction (LDR)



LDR conditions: 1X Ligase Buffer (20 mM Tris (pH 7.5), 20 mM KCl, 1 mM DTT, 10 mM MgCl₂, 0.1% IGEPAQ®, 0.01 mM Cofactors, 0.50 μM Donor probe, 0.25 μM Acceptor probe, ssDNA template (50 nM), 4 Units *Pfu* DNA Ligase, 20 μL. Cycling conditions: 94°C (2 min), Cycle 20x [94°C (30 sec), 65°C (4 min)].

Conclusion

- 1) Chemical modifications made to the ATP cofactor, donor probes and acceptor probes provide improved ligation specificity for T4 DNA ligase.
- 2) Optimal chemical modifications to acceptor probes vary depending on the DNA ligase family (ATP or NAD⁺).
- 3) Cofactor modifications provide the greatest improvement to ligation fidelity of T4 DNA ligase.
- 4) The lead ATP cofactor demonstrated the most consistent increase in ΔCq's in real-time PCR.
- 5) LDR reactions employing chemically modified ligation components yield more specific match vs. mismatch discrimination.

Acknowledgements:

Terry Beck, Angela Tenenini, Anton McCaffrey, Julie Powers, Elizabeth Hill


NIH Grants:

1R43GM085860-01, 2R44GM085860-02 & 5R44GM085860-03

For further information, please contact Natasha Paul, npaul@trilinkbiotech.com

Figure 6


Comparison of acceptor probe modifications using ATP or NAD⁺ ligases



Ligation conditions: 1X Ligase Buffer (50 mM Tris (pH 7.5), 10 mM DTT, 10 mM MgCl₂, 1 mM ATP), Ligation Probes (1 μM), ssDNA Template (1 μM), T4 DNA Ligase or 2 Units *E.Coli* DNA Ligase, 20 μL. T4 DNA Ligase Thermal cycling conditions: 95°C (3 min), 4°C (3 min), add Ligase, 22°C (20 min), 65°C (10 min). *E.Coli* DNA Ligase Thermal cycling conditions: 95°C (3 min), 4°C (3 min), add Ligase, 16°C (20 min), 65°C (20 min).

Figure 7

Investigation of chemically modified donor probe modifications



Ligation conditions: 1X Ligase Buffer (50 mM Tris (pH 7.5), 10 mM DTT, 10 mM MgCl₂), Ligation Probes (1 μM), ssDNA Template (1 μM), T4 DNA Ligase 2 Units, 20 μL. Thermal cycling conditions: 95°C (3 min), 4°C (3 min), add Ligase, 16°C (20 min), 65°C (10 min).