

# METABOLISM OF EIGHT MODEL PHARMACEUTICAL COMPOUNDS IN RAT- and HUMAN- HEPATOPAC® VERSUS LIVER MICROSOMES AND SUSPENSION HEPATOCYTE PLATFORMS



Julius O. Enoru<sup>1</sup>, William DeMaio<sup>1</sup>, Adiba Watanyar<sup>1</sup>, Kenneth Draper<sup>1</sup>, Amanda Moore<sup>2</sup> & Okey Ukairo<sup>2</sup>

1. Ricerca Biosciences LLC, 7528 Auburn Rd, Concord, OH 44077  
2. Hepregen Corporation, 200 Boston Ave, Medford, MA 02155



## ABSTRACT

Metabolism evaluations for compounds in early development are typically conducted *in vitro* using liver microsomes and suspension hepatocytes. While these systems have served the drug development community for many years, they have limitations and result in infrequent “metabolite surprises”, a term used to describe the observation of human metabolites that were previously undetected in preclinical studies. Such surprises can have a significant impact on the timeline of a drug development program and has triggered investigations to develop *in vitro* systems that more accurately predict human *in vivo* metabolism.

In this study, the metabolite profiles of eight model pharmaceutical compounds with various biotransformation reactions were investigated in parallel using a functionally stable model of primary hepatocytes [micropatterned co-cultures (MPCCs), referred to as Hepatopac™] and the traditional liver microsomal and suspension hepatocyte platforms (Khetani & Bhatia, 2006). Incubations were performed for 0 and 60 min in liver microsomes; 0 and 2 hr in suspension hepatocytes; and 0, 4 hr, 2 days and 7 days in Hepatopac™. Metabolites were identified by LC/MS employing a narrow range of chromatographic conditions, which were representative of drug metabolism screening parameters used in an early drug development setting.

Our results show that Hepatopac™ are metabolically more active than the traditional platforms, as judged from levels of test article disappearance and formation of metabolites. Formation of metabolites from non-P450 mediated reactions, especially UGTs, and secondary metabolites appear to be more predominant in Hepatopac™.

Results from this investigation also suggest that use of Hepatopac™ for metabolite profiling studies will be broadly useful in minimizing “metabolite surprises”. Use of Hepatopac™ platform also provides cells that maintain metabolic activity over an extended period in culture; a difference that likely explains the increased abundance of secondary metabolites in this system.

## INTRODUCTION

Regulatory guidance specifies safety testing of new human metabolites which represent greater than 10% of the parent drug-related components in circulation. Early *in vitro* metabolite profiling studies can provide a clue regarding any potential differences in metabolite profiles between humans and animals. Such *in vitro* evaluations are typically conducted using human liver microsomes and suspension hepatocytes. These systems however have their limitations, sometimes resulting in “metabolite surprises” - where metabolites which were not predicted or observed *in vitro* or *in vivo* during pre-clinical studies, are observed in the clinic. Such surprises can have a significant impact on the timeline of a drug development program and has prompted investigations aimed at developing new *in vitro* systems and models that can more accurately predict human *in vivo* metabolism.

This poster presents our comparative metabolite profiling analysis of eight model pharmaceutical compounds with various biotransformation reactions using a functionally stable model of primary hepatocytes [micropatterned co-cultures (MPCCs)] in parallel with the traditional liver microsomal and suspension hepatocyte systems.

### Current In Vitro Platforms and Limitations


- Liver microsomes fortified with cofactors and cytosol- limited non-P450 enzymes
- Suspension hepatocytes - Loses phenotype in culture
- Liver S9 fortified with cofactors - Limited amount of phase II enzymes
- Recombinant enzymes - “More artificial” than other *in vitro* systems

Note: All systems show a decrease in P450 enzyme activity with time in culture

### Micropatterned Platform


- Micropatterned co-cultures of hepatocytes and other cells
- 3-D context
- Cell-cell interactions
- Cell-Matrix interactions
- Maintains metabolic activity over extended periods in culture

**Figure 1: Zaleplon – Metabolite Profiles in Human LM, Hepatocytes and Micropatterned Platform**




HepatoPac™ and cytosol fortified LM predicted actual aldehyde oxidase (AO) species difference between humans and rats as expected

**Figure 2: Zaleplon – Metabolite Profiles in Human Plasma 1 hr after 20 mg dose**  
(Reference: DeMaio et al., 1994.)




HepatoPac consistent with human plasma radiochromatogram, except for 278, which was present as an additional metabolite. 5-oxo ZAL, a minor metabolite also observed in rat bile and urine, was not observed with rat HepatoPac. Desethyl ZAL metabolite (278) was also the major metabolite in dog plasma.

**Figure 3: Lorcaserin – Metabolite Profiles in Human LM, Hepatocytes and Micropatterned Platform**




A major human metabolite carbamoyl glucuronide (CG) is generated by HepatoPac™ but not predicted by LM and hepatocytes.

**Figure 4: Ziprasidone – Metabolite Profiles in Human LM, Hepatocytes and Micropatterned Platform**




HepatoPac™ and hepatocytes confirmed aldehyde oxidase and S-methyl transferase metabolism that was not observed in Liver MIC+CYT

**Figure 5: Buspirone – Metabolite Profiles in Human LM, Hepatocytes and Micropatterned Platform**



All systems generated numerous oxidation products. HepatoPac™ metabolism was more extensive than other systems.

**Figure 6: Nefazodone – Metabolite Profiles in Human LM, Hepatocytes and Micropatterned Platform**




Human HepatoPac™ closest to human plasma profile than the traditional in vitro platforms

This study was jointly sponsored by Ricerca Biosciences, LLC, Concord, OH, & Hepregen Corporation, Medford, MA


ISSX 2014

**Figure 7: XK-469 – Metabolite Profiles in Human LM, Hepatocytes and Micropatterned Platform**




XK-469 was minimally metabolised in all systems, consistent with human plasma profile

**Figure 8: Carbazeran – Metabolite Profiles in Human LM, Hepatocytes and Micropatterned Platform**



Human HepatoPac™ showed more glucuronidation than the other *in vitro* systems

**Figure 9: Zomepirac – Metabolite Profiles in Human LM, Hepatocytes and Micropatterned Platform**



Minimal metabolism observed in all platforms

## MATERIALS and METHODS

### Model Compounds

- Zaleplon (ZAL)
- Ziprasidone (ZIP)
- Buspirone (BUS)
- Nefazodone (NEF)
- Lorcaserin (LOR)
- XK-469
- Carbazeran (CAR)
- Zomepirac (ZOM)

### Experimental Design:

- Liver Microsomes (LM) with Cyt - 1 hr incubation in presence of NADPH and UDPGA
- Suspension Hepatocytes - Incubation time: 2 hr
- HepatoPac™ - Incubation time: 4 hr, 2 days, 7 days
- Compound concentration = 10 μM
- Note: Suspension hepatocytes and HepatoPac™ were derived from the same lot of primary hepatocytes.

### LC/MS Analysis:

Use same HPLC column and mobile phase for all compounds  
Determine DP (Q1MS scan) and CE (MS/MS) for each parent drug  
Determine mobile phase gradient to get parent drug retention time in 10 to 15 min range, starting gradient at 5% organic.  
Phenomenex Luna C18(2) 100 x 2.00 mm, 3 micron  
Mobile Phase A = 10 mM NH<sub>4</sub> acetate in water, pH 4.5  
Mobile Phase B = Acetonitrile  
Run time = 30 min  
LC/MS (+) ESI mode AB Sciex API4000, Q1 scan, undiluted and 10x-dilution  
LC/MS/MS (+)ESI mode, product ion scans for selected metabolites

## SUMMARY OF RESULTS

- Zaleplon - Actual AO species difference between humans and rats predicted by HepatoPac™ and cytosol fortified liver microsomes as expected
- Lorcaserin - A major human metabolite (carbamoyl glucuronide) is generated by HepatoPac™ that was not predicted by LM or HEP
- Ziprasidone - Aldehyde oxidase and S-methyl transferase metabolism confirmed in hepatocytes and HepatoPac™ that was not observed in Liver MIC+CYT
- Buspirone - Numerous oxidation products observed in all systems as expected. HepatoPac™ metabolism more extensive than other systems
- Nefazodone - Human HepatoPac™ closest to human plasma profile than the traditional *in vitro* platforms
- XK-469 - Minimal metabolism in all systems consistent with minimal metabolism indicated by human plasma profile
- Carbazeran - More glucuronidation observed in human HepatoPac™

## CONCLUSION

- HepatoPac™ maintains enzyme activity for longer periods.
- HepatoPac™ provides more extensive metabolism than other platforms and therefore may serve as a useful system for predictive toxicology evaluations.
- Metabolic stability in HepatoPac™ appears more predictive of actual human *in vivo* metabolic stability than obtainable with suspension hepatocytes.
- UGT activity higher with HepatoPac™ than traditional platforms
  - N-carbamoyl glucuronidation activity high
  - Acyl glucuronidation activity low.
- Aldehyde oxidase (AO) species differences readily apparent (ZAL and CAR)
- Metabolite profiling with HepatoPac™ is a service that would be valuable to our drug development clients.

## REFERENCES

1. Khetani S.R. and Bhatia S.N. Engineering Tissues for In Vitro Applications. Current Opinion in Biotechnology, Vol 17, p524-531 (2006).
2. Z. Tong, A. Chandrasekaran, W. DeMaio, R. Espina, W. Lu, R. Jordan, and J. Scatina. 2010. Metabolism of Varicaserin in Mice, Rats, Dogs, Monkeys and Humans. *Drug Metab Dispos*, 38 (12), 2266-2277.
3. W. DeMaio, G. Morton, I. Chaudhary, R.G. Kelly. 1994. Identification of metabolites of a novel sedative/hypnotic agent by FAB-MS/MS and HPLC/ionspray-MS of radio and stable isotope labeled samples. Proc. 42nd ASMS Conference, Chicago, IL.

HepatoPac™ is a trademark of Hepregen Corporation