



qNMR is the “gold standard” for compound quantification and has wide application. The quantitative analysis of mixtures poses special difficulties, but is a very worthwhile objective.

We describe a simple, robust approach to mixtures analysis where each component can show a distinct NMR signal, or one related to another component. The method is applicable to a wide range of foods, cosmetics, body fluids, *etc.* This is called *Simple Mixtures Analysis (SMA)*.

## Basic processing: GSD




This trial NMR spectrum (red) is challenging for deconvolution  
GSD still deconvolutes the peaks (blue) and determines areas



The two doublet-doublets of Citric Acid here overlap heavily with other mixture components in wine

Multiplet areas differ only by <0.5% using GSD

## SMA examples



**Results:**  
Compound Result (mg/L)  
004 Acetic acid (2.08 ppm, s) 204.00  
003 Salicylic acid (7.88 ppm, dd, 7.9 Hz, 1.9 Hz) 259.40  
002 Ascorbic acid (4.95 ppm, 1.9 Hz) 475.04  
001 Acetylsalicylic acid (7.97 ppm, dd, 7.8 Hz, 1.7 Hz) 907.71  
005 Citric acid (2.98 ppm, d, 15.7 Hz) 3823.47

## SMA - Overview


### Library experiments

A *protocol*, or *Method* to measure component concentrations in a mixture must first be developed and validated by the user. Data from 1D or 2D experiments can be used.

Application of this analysis to routine samples is then significantly expedited, using a work-flow that combines full automatic software analysis and user checking.

*Integration data is extracted from multiplet peaks in each new spectrum, and a user-specified equation determines the concentration or weight%.*

The Experiment was redesigned now to individually assess the concentration based on all component multiplets in Aspirin Plus C



### The workflow improves speed and quality

#### Error checking

Multiplet type,  $J_s$ , and integral ratios can be specified. If outside acceptable ranges, a warning is produced


#### Manual refinement

Incorrectly selected multiplets can be manually adjusted for a correct analytical result, and the concentrations corrected

#### Reporting

A simple output summary can be pasted on the spectrum, or an XML file written for customized output or DB inclusion.

## Method development: Experiments




An *Experiment* is made up of a definitions for a sample reference, and known components:

- PPM ranges and specified multiplets define each species by its area
- If target multiplets are not correctly found then an error message can be shown


The *equation editor* flexibly specifies calculation of a concentration

- A mixture reference can be used to derive a  $CCF_i$  which can then be used to determine other component concentrations
- Spectral math: use areas and sample-specific data from >1 experiment
- 1D and 2D datasets
- Integral data from the spectra

### SMA of amino acids in serum



## Component peaks



### Peak Pattern Recognition

Component multiplets can be recognized using conventional methods, or with a novel pattern recognition algorithm.

A Reference peak pattern is found, and searched for in the mixtures spectrum. A scoring system is used, which also allows inexact but correct matches to be found.

We have designed a versatile, functional program for the quantitation of simple mixture components by NMR. The user can use a flexible equation editor to have access to integrations from one or more 1D or 2D spectra, and thereby design an analysis. Multiplets can be automatically checked and the analysis adjusted manually. Reporting can be fully customized.

The software is applicable to a wide range of real-world mixtures such as foods, nutraceuticals, *etc.*

### Acknowledgements

The authors thank Dr Silvia Mari (Research4rent) for the serum spectrum, and Dr Susanne Klein (Institut Heidge) for the Aspirin Plus C spectra and helpful discussions.