

# AUTOMATED TYROSINE KINASE INHIBITOR CARDIOTOXICITY ASSAY IN ZEBRAFISH

Olaia Holgado<sup>1</sup>, Juan Maria Virto<sup>1</sup>, Izaskun Ibarbia<sup>1</sup>, Patrice Dubreuil<sup>2</sup>, Didier Pez<sup>2</sup>, Ainhua Letamendia<sup>1</sup>, Martine Humbert<sup>2</sup>, Alain Moussy<sup>2</sup>, Carles Callol-Massot<sup>1</sup>.

<sup>1</sup>BIOBIDE S.L., Paseo Mikeletegi 58, Donostia, 20009, Spain ([WWW.BIOBIDE.COM](http://WWW.BIOBIDE.COM))

<sup>2</sup>AB Science, France ([www.ab-science.com](http://www.ab-science.com))

## Abstract

The zebrafish embryos have recently gained relevance in biomedical research thanks to some of its characteristics including embryo transparency, small size, ease of manipulation and possibility to evaluate different internal organs avoiding invasive methodologies. Combined with the possibility to adapt the model with an automatic device and the reduced cost associated to each assay, the model is an ideal killer experiment in early phases of drug discovery as well as a novel method to increase the selection arguments to reduce the candidates to enter into the Drug Development processes.

Cardiotoxicity is one of the most important reasons for drug attrition during the process of Drug Development. Evaluation of cardiotoxicity and especially HERG channel inhibition is described in regulatory guidelines, but limitations demands the development of new complementary assays that can also evaluate the heart function from a holistic point of view. Biobide has set up a novel in vivo automated platform that allows testing compounds in zebrafish embryos.

To evaluate and validate the quality of the analysis system, the model and the value of the information, we have used a panel of blind-coded Tyrosine kinase inhibitors that had been previously described in other in vitro and in vivo assays. The results indicate that our automated method provides with high informative and complementary data that can significantly improve the process of selection of new candidates with low or no cardiotoxicity.

## Method

Embryos expressing a green fluorescent protein under a heart specific promoter were collected at 0 hours post fertilization (hpf) and raised with E3 media at 28.5C in an incubator. Once they reached 48 hpf, embryos were poured into an embryo sorter and sorted based on the expression of fluorescence and size of the embryo. Embryos were incubated with each compound at 3 different concentrations 5, 10, 50  $\mu$ M and analyzed at 3 and 24 hours of incubation. All the processes were done automatically in the platform of Biobide (figure 1). Heart beat alterations were automatically analyzed with algorithms specially developed (figure 2) and a final comparison with other methods was done with data provided by AB Sciences.




Figure 1. Images obtained with the automated system. Heart beat analysis is performed using videos of fluorescent transgenic hearts of 48 hpf zebrafish embryos.




Figure 2. Automated workflow in a Validated Screening Platform Once zebrafish embryos are obtained from massive crossing, they are dispensed in 96 wells plate automatically.

After this, compounds are added at the specified concentration and a robotic arm takes the plate to the image analysis station. In house developed softwares allow choosing the safest /most efficient compounds.

## Results

### Description of the Tyrosine Kinases

Tested drug	Clinic General Description	Cardiac effects	Incidence cardiac event in clinic
SORAFENIB	Is an antitumoral agent. Among the adverse effects of Sorafenib cardiovascular effects can occur, such as cardiomyopathy and hypertension to QT prolongation. Sorafenib also causes disthyroidism that aggravates cardiotoxicity. The registratroy trials of Sorafenib reported cardiotoxicity in approximately 10% of the patients, but more recent observational studies of the general population show that cardiac events may occur in as many as 40% of patients.	Hypertension, Hypothyroidism that aggravates cardiotoxicity, mytochondrial dysfunction, ischemia	17% Hypertension, Warning Full prescribing information
SUNITINIB	Sunitinib is a multitargeted TKI used to extend survival of patients with metastatic renal-cell carcinoma and gastrointestinal stromal tumours, but concerns have arisen about its cardiac safety, with reported rates of cardiotoxicity are as high as 33.8% .	QT prolongation, Cardiomyopathie and Hypertension, mytochondrial dysfunction, ischemia	18% patients, 24% Hypertension, Warning Full prescribing information
DASATINIB	Dasatinib is a Src family TKI used in patients with chronic myelogenous leukemia (CML). There is limited data on Dasatinib and cardiotoxicity. However, QT prolongation has been observed in Phase II trials of Dasatinib, and severe cardiac toxicity has occasionally been reported with Imatinib.	QT prolongation	Low to moderate
IMATINIB	Imatinib is a TKI used to treat certain types of cancer. The most common side effects include weight gain, reduced number of blood cells (neutropenia, thrombocytopenia, anemia), headache, edema, nausea, rash, and musculoskeletal pain. Severe cardiotoxicities are occasional.		No or low rate
NILOTINIB	Nilotinib is used in treatment of leukemia. Nilotinib prolongs the QT interval, and sudden deaths have been reported rarely.	QT prolongation	Warning FPI
PAZOPANIB	Pazopanib is a second-generation multitargeted TKI that has exhibited antiangiogenic and antitumor activity. Pazopanib has been shown to prolong the QTC interval (2% of the cases in clinical studies). The concurrent use of Pazopanib with other agents that prolong the QTC interval may result in potentially life-threatening cardiac arrhythmias, including Torsades de Pointes.	QT prolongation and Hypertension	37% Hypertension, Warning Full prescribing information
MASITINIB	Masitinib is a selective inhibitor of c-Kit and PDGFR registered in Europe (Masivet®) and the USA (Kinavet®) for the treatment of canine mast cell tumor, a common cutaneous cancer in dogs. Masitinib is being evaluated in nine phase III studies in human medicine, including pancreatic cancer, GIST, metastatic melanoma, multiple myeloma, mastocytosis, severe persistent asthma, rheumatoid arthritis, Alzheimer's disease, and multiple sclerosis		-
BOSUTINIB	Bosutinib functions as a dual inhibitor of Src and Abl kinases, with little if any inhibitory effect over PDGF-R and Kit. It is used in the treatment of leukemia.		-
PONATINIB	Ponatinib is a multi-targeted TKI, being its primary target BCR-ABL and it is an experimental oral drug candidate for the treatment of chronic myeloid-leukemia.		-

Cardiotoxicity classification is based on the use of the euclidian distance algorithm.




Figure 3. Description of the drugs tested in the zebrafish cardiotoxicity platform in Biobide, related to their cardiotoxic known effects in humans.

Tyrosine kinases potential cardiotoxicity was assessed in zebrafish embryos automatically in the platform at 3 and 24 hours post incubation

Drugs	Solubility	Concentration (μM)						Final conclusion
		3h treatment			24h treatment			
	Water μM	5	10	50	5	10	50	
SORAFENIB	8	Br	Br	Br	D	D	D	Bradycardia at 3h, acute toxicity at 24h
SUNITINIB	150	NE	NE	NE	NE	Br/A <sup>1</sup>	Br/A <sup>1</sup>	<sup>1</sup> Dysfunction Ventricle and Bradycardia
								zERG blocker
DASATINIB	200	NE	NE	NE	NE	NE	Br	Moderate Bradycardia
IMATINIB	200	NE	NE	NE	NE	NE	NE	No cardiotoxicity
NILOTINIB	2	NE	NE	Br	NE	NE	NE*	Moderate Bradycardia
PAZOPANIB	31	NE	NE	NE	NE	NE	NE	No cardiotoxicity
MASITINIB	100	NE	NE	NE	NE	NE	NE	No cardiotoxicity
BOSUTINIB	16	NE	NE	D	NE	NE	D	No clear cardiotoxicity but acute toxicity and heart failure.
PONATINIB	15	NE	NE	NE	D/A <sup>1</sup>	D/A <sup>1</sup>	D/A*	<sup>1</sup> Dysfunction Ventricle
								zERG blocker

Figure 5. Results of the tested drugs in the in vivo zebrafish automated platform in Biobide (toxicities in red)

No effect (NE); Bradycardia (BR); Arrhythmia (A) ; Death (D); \* Precipitation observed

The zebrafish assay is a good tool to predict cardiotoxicity in humans and can complement in vitro assays.

Drugs	hERG# Inhibition (%)	hCM Prolif (IC50 $\mu$ M)	hROS prod*	Mito Resp <sup>2</sup> (%)	Incidence of cardiac event in clinic	Zebrafish
SORAFENIB	9	>10	-21	33	17% HTN Warning FPI	Detected
SUNITINIB	52 (0.27 $\mu$ M)*	9.5	331	93	18% patients 24% HTN Warning FPI	Detected
DASATINIB	23 (14 $\mu$ M)*	<1	167	134	Low to mod rate	Detected
IMATINIB	15	>10	15	160	No or Low rate	Not detected
NILOTINIB	40 (0.42 $\mu$ M)*	3	240	35	Warning FPI	Detected
PAZOPANIB	33	>10	39	ND	37% HTN Warning FPI	Not detected
MASITINIB	23	>10	11	100	-	Not detected
BOSUTINIB	96	>10	7	56	-	Not classified because of dose
PONATINIB	86	<1	ND	125	-	Detected

Figure 6. List of TKIs and potential cardiac activity in vitro. (cardiotoxicites in red) # % inhibition at 3 $\mu$ M Predictor hERG (Invitrogen) ;\*Patch clamp IC<sub>50</sub>  $\mu$ M; human cardiomyocytes hCM ; hROS ROS production by treated cardiomyocytes ; \*\* % Respiration expressed as a ratio to a non toxic compound in isolated murin hearts

## Conclusions

- This zebrafish test classified correctly 7 out of the 9 compounds that were tested.
- All compounds that were reported to be cardiotoxic in clinic were also positive in the zebrafish assay with the exception of Pazopanib that was not identified in any of the assays. The main cardiotoxic activity is based on hypertension, an endpoint that was not evaluated in the zebrafish assay and that could be implemented in future assays.
- Testing of compounds at 3 and 24 hours post incubation was informative and complementary. The incubation at 24 hours allowed the evaluation of compounds inducing alterations in metabolism or cardiomyocyte survival such as with the Sunitinib.
- Testing compounds in zebrafish embryos can help establishing the therapeutic window when combining different concentrations and specific analysis with general toxicity analysis.
- Bosutinib induced embryo toxicity observed as heart failure and reported as embryo dead. No cardiotoxicity was reported at lower doses in accordance with the results observed in clinic. New analysis at concentrations between 10 and 50  $\mu$ M would be required to finally classify this compound.
- The zebrafish cardiotoxic data is complementary to other already developed assays to evaluate cardiotoxicity.

## Bibliography

Langheinrich U, Vacun G, Wagner T. Zebrafish embryos express an orthologue of HERG and are sensitive toward a range of QT-prolonging drugs inducing severe arrhythmia. Toxicol Appl Pharmacol. 2003 Dec 15;193(3):370-82.

Milan DJ, Peterson TA, Ruskin JN, Peterson RT, MacRae CA. Drugs that induce repolarization abnormalities cause bradycardia in zebrafish. Circulation. 2003 Mar 18;107(10):1355-8.

Burns CG, Milan DJ, Grande EJ, Rottbauer W, MacRae CA, Fishman MC. High-throughput assay for small molecules that modulate zebrafish embryonic heart rate. Nat Chem Biol. 2005 Oct;1(5):263-4. Epub 2005 Sep 18.