

# The advancement of feeder-independent culture systems for undifferentiated human pluripotent stem cells towards greater regulatory compliance

Jennifer Moody<sup>1</sup>, Jennifer Antonchuk<sup>1</sup>, Debbie King<sup>1</sup>, Andrea Eskander Afshari<sup>1</sup>, Heather Drew<sup>1</sup>, Alexandra Blak<sup>1</sup>, Clive Glover<sup>1</sup>, Cindy Miller<sup>1</sup>, Terry E. Thomas<sup>1</sup>, Allen C. Eaves<sup>1,2</sup> and Sharon A. Louis<sup>1</sup>

<sup>1</sup> STEMCELL Technologies Inc., Vancouver, BC, Canada

<sup>2</sup> Terry Fox Laboratory, BC Cancer Agency, Vancouver, BC, Canada

## Introduction

Defined and feeder-independent cell culture systems provide a platform for greater reproducibility and standardization in human pluripotent stem cell (hPSC) research. mTeSR<sup>®</sup>1 has become the most widely cited feeder independent system for the culture of undifferentiated human ES and iPSC (Table 1). As the field develops potential therapeutic applications for hPSC-derived cells, it is increasingly important that media products are manufactured to meet applicable regulatory compliance standards. Therefore, we have developed two new products for the expansion of undifferentiated hPSCs; mTeSR<sup>®</sup>1 manufactured in GMP facility (mTeSR<sup>®</sup>1-GMP) and animal protein-free TeSR<sup>™</sup>2. Here we show data to support their use for maintenance of high quality undifferentiated hPSC cultures.

Figure 1: The TeSR<sup>™</sup> System




Table 1. Ten select publications using mTeSR<sup>®</sup>1 in 2009 for a variety of applications

### Derivation and/or Culture of iPSC

Chan et al., 2009 Nat Biotechnol. Nov;27(11):1033-7. Live cell imaging distinguishes bona fide human iPS cells from partially reprogrammed cells.

Eminli et al., 2009 Nat Genet 4:968-976. Differentiation stage determines potential of hematopoietic cells for reprogramming into induced pluripotent stem cells.

Sun et al., 2009 PNAS Sep 15;106(37):15720-5. Feeder-Free derivation of induced pluripotent stem cells from adult human adipose stem cells.

### Derivation of hESC

Eremeev et al., 2009 Dokl Biol Sci 426:293-295. Derivation of a novel human embryonic stem cell line under serum-free and feeder-free conditions.

### Cell Source for Differentiation of hPSC

Yu et al., 2009 PLoS ONE 4(9): e7040. nAChRs mediate human embryonic stem cell-derived endothelial cells: proliferation, apoptosis, and angiogenesis.

Lu et al., 2009 Exp Hematol 37:924-936. Enhanced generation of hematopoietic cells from human hepatocarcinoma cell-stimulated human embryonic and induced pluripotent stem cells.

### hESC Characterization

Zeng et al., 2009 Stem Cells Oct;27(10):2435-45. Lack of ABCG2 expression and side population properties in human pluripotent stem cells.

Kolle et al., 2009 Stem Cells Oct;27(10):2446-56. Identification of human embryonic stem cell surface markers by combined membrane-polysome translation state array analysis and immunotranscriptional profiling.


### Alternate Culture Systems

Hakala et al., 2009 Tissue Eng Part A. Jul;15(7):1775-85. Comparison of biomaterials and extracellular matrices as a culture platform for multiple, independently derived human embryonic stem cell lines.

Oh et al., 2009 Stem Cell Research Mar 4. [Epub ahead of print]. Long term microcarrier suspension cultures of human embryonic stem cells.


## Results

Figure 2. Expansion of undifferentiated cells in mTeSR<sup>®</sup>1-GMP is comparable to their expansion in mTeSR<sup>®</sup>1




a) Cells were cultured in mTeSR<sup>®</sup>1 and mTeSR<sup>®</sup>1-GMP on Matrigel<sup>™</sup> for 5 passages with minimal selection of differentiated colonies (<10% per passage). Expansion was measured by enumerating clumps at harvest/clumps seeded (n=2, data pooled from 4 passages each  $\pm$  SEM). b) Oct3/4 expression measured by FACS was >95% at the end of 5 passages for H9 cells cultured in mTeSR<sup>®</sup>1 (orange, solid line) and mTeSR<sup>®</sup>1-GMP (pale orange, dotted line). Secondary antibody only control for mTeSR<sup>®</sup>1-GMP = grey, solid line.

Figure 3. Long-term culture of undifferentiated cells in animal protein-free TeSR<sup>™</sup>2 affords robust expansion and maintenance of pluripotency marker expression




a) H1 and H9 were maintained in TeSR<sup>™</sup>2 on Matrigel<sup>™</sup> for > 25 passages with minimal selection of differentiated colonies (<10% per passage). Average ( $\pm$  1 SEM) expansion per passage was measured as described in Figure 1 (24 passages for H9; 22 passages for H1). b) FACS analysis of Oct3/4, SSEA-4 and Tra-1-81 after >25 passages in TeSR<sup>™</sup>2.

Figure 4. hPSC maintained in TeSR<sup>™</sup>2 are functionally pluripotent




H9 cells isolated after 11 passages in TeSR<sup>™</sup>2 on Matrigel<sup>™</sup>, and then injected into mice, generated teratomas over a period of 9 weeks. a, b) Histological analysis of H&E stained sections of the resulting teratomas revealed tissues representative of all three germ layers as indicated. gut epithelia = endoderm; neural=ectoderm; cartilage, muscle=mesoderm

Figure 5. hPSC cultured in TeSR<sup>™</sup>2 maintain a normal karyotype with prolonged passaging



G band karyotype analysis of a) H9 cells (normal, XX) after 22 passages in TeSR<sup>™</sup>2 on Matrigel<sup>™</sup> and b) H1 cells (normal, XY) after 19 passages in TeSR<sup>™</sup>2 on Matrigel<sup>™</sup>.

Figure 6. Undifferentiated growth of hPSC in TeSR<sup>™</sup>2 on a defined matrix



a) H9 and b) H1 cells were cultured in TeSR<sup>™</sup>2 on recombinant human vitronectin (R&D Systems) for 5 passages using mechanical dissociation of clumps without enzyme at each passage. Cells cultured in these conditions form tightly packed colonies and maintain their characteristic undifferentiated morphology. c) Average ( $\pm$  1 SEM) expansion per passage was measured as described for Figure 1 (4 passages for each line). FACS analysis at the end of 5 passages reveals >95% SSEA4 (d) and >97% SSEA3 (e) expression by both H1 and H9 cells. Secondary antibody only control for H1 = gray, solid line, H1=pale orange, dotted line, H9=orange, solid line.

## Conclusions

- mTeSR<sup>®</sup>1 manufactured in a GMP facility provides equivalent performance to mTeSR<sup>®</sup>1.
- TeSR<sup>™</sup>2 is an animal protein-free formulation that can be used with Matrigel<sup>™</sup> or more defined matrices to maintain pluripotent and karyotypically normal hPSC.
- These two new products represent STEMCELL's ongoing commitment to provide ongoing regulatory compliant reagents to facilitate the translation of basic research to pre-clinical studies.



THE CELL EXPERTS<sup>™</sup> | WWW.STEMCELL.COM

IN NORTH AMERICA TOLL-FREE T. 1 800 667 0322 TOLL-FREE F. 1 800 567 2899  
T. 1 604 877 0713 F. 1 604 877 0704 E. INFO@STEMCELL.COM  
TOLL-FREE T. 00 800 7836 2355 TOLL-FREE F. 00 800 7836 2300  
T. +33 (0)4 76 04 75 30 F. +33 (0)4 76 18 99 63 E. INFO.EU@STEMCELL.COM  
IN AUSTRALIA TOLL-FREE T. 1 800 060 350 F. +61 (03) 9338 4320 E. INFO.AUS@STEMCELL.COM  
IN SINGAPORE T. (65) 6776 7754 F. (65) 6776 7114 E. INFO.SG@STEMCELL.COM