

Sensitive Identification of Phosphopeptides in Brain Tissue using 2D-NanoLC-ESI-MSⁿ

Jenny Samskog, Henrik Wadensten, and John Flensburg
GE Healthcare/ Amersham Biosciences AB, Uppsala, Sweden

Introduction

One of the most important post-translational modifications is phosphorylation of serine, threonine or tyrosine residues. Detection of phosphorylation sites by mass spectrometry in proteins extracted from biological material is complicated by low abundance, low stoichiometry, and poor ionization of phosphopeptides [1]. In this work, a biocompatible nano liquid chromatography (LC) system, Ettan™ MDLC, was used for separating tryptic peptides from brain tissue by cation exchange (SCX) to enrich the phosphopeptides followed by reversed-phase chromatography (RPC). The phosphopeptides were detected by neutral loss MS.

Methods

Mouse brain tissue was trypsin digested and analysed using Ettan MDLC (Amersham Biosciences) coupled to a Finnigan LTQ linear ion trap (Thermo ElectronCorp.). 40 µg of sample was injected onto a 2.1 * 250 mm SCX column (BioBasic, Thermo Electron) and eluted with a linear salt gradient (A: 20 mM Citric acid, 25% CH₃CN, B: A + 1 M NH₄Cl) where fractions were collected (Fig.1).

The fractions were injected onto a trap column (Zorbax, Agilent, 0.3*5 mm) and RPC separation was performed on a 0.075 * 150 mm Zorbax column (Agilent). Two sets of trap/separation columns




Fig 1. The MDLC off-line configuration.

The MS method consisted of a cycle combining one full MS scan with three MS² events (25% collision energy) followed by a MS³ event (35% collision energy) that was triggered upon detection of -98, -49, or -32.7 Daltons from the precursor (neutral loss of phosphoric acid, charge states 1+, 2+, and 3+). Dynamic exclusion duration was set to 30 s. The MS² and/or MS³ spectra from all the runs were searched using TurboSEQUEST protein identification software (Thermo Electron). Modifications were set to allow for the detection of oxidized Met (+16), carboxyamidomethylated Cys (+57), phosphorylated Ser, Thr, and Tyr (+80), and dehydrated Met and Thr (-18). The peptide matches were filtered based on cross-correlation scores (Xcorr) of 1.5, 2, and 2.5 for charge states 1+, 2+, and 3+, respectively.

1. Ficarro, S. B. *et al.* Phosphoproteome analysis by mass spectrometry and its application to *Saccharomyces cerevisiae*. *Nat. Biotechnol.* **20**, 301–305 (2002).
2. Beausoleil, S. A. *et al.* Large-scale characterization of HeLa cell nuclear phosphoproteins. *Proc. Natl. Acad. Sci. USA* **101**, 12130–12135 (2004).

Results

The neutral loss MS method was optimised for single dimension RPC separation. The detection limit was then shown to be less than 1 fmol for a phosphopeptide in a protein digest standard. In this work another separation dimension was added to the system to both increase the chromatographic resolution in the system and to concentrate the phosphopeptides [2] by SCX.




Fig 2. UV trace from SCX separation. The green bars indicate the relative amount of phosphopeptides that were identified in the fractions.

The SCX separation of the tryptic digest can be seen in Fig. 2. The phosphopeptides eluted between 8-17% NH₄Cl and were limited to 11 fractions of totally 31 analysed. Fraction, b7, contained most of the identified phosphopeptides, namely 15 peptides. Most of the phosphopeptides were only found in one fraction which indicates that the size of the fractions (collected every 30 second) correlated well with the peak width.




Fig 3. Base peak ion chromatogram from fraction b10 and all MS³ events.

The TurboSEQUEST searches resulted in some false positives, often due to incorrectly assigned charge states for peptides eluting late in the RPC run, see Fig. 3. Some examples of peptides that were found to be phosphorylated in both MS² and MS³ database searches are shown in Fig 4. The neutral loss ion is apparent in the MS² spectrum and the peptide is sequenced from the MS³ spectrum. Some of the identified phosphopeptides are shown in Table 1.




Fig 4. MS² and MS³ spectra from ApSGQAFELILSPR (from Stathmin 1)

The developed strategy for confident analysis of phosphopeptides in complex mixtures is summarised below:

1. 2D LC (SCX/RPC)
2. MS³ on all peptides that loses phosphoric acid (neutral loss)
3. TurboSequest searches on all MS³ spectra (-18@ST)
4. Manual confirmation of charge state and that neutral loss dominates MS/MS spectra.
5. Further confirmation by MS² searches of +80@STY




Fig 5. A tyrosine-phosphorylated peptide originating from seryl-aminoacyl-tRNA synthetase 2, Xcorr: 2.6

Table 1. Some of the phosphopeptides that were identified in both MS² and MS³ using TurboSEQUEST.

Protein	Sequences ^a	MH+	Xcorr ^b	SCX fraction
Microtubule associated protein 1B	pPEEGGVVESEK	1570	4.3	b6
	RpGKETETK	1386	3.4	c2, c3
	ADpSRESKPAIK	1284	3.2	c2
Similar to Microtubule-associated protein 2	CYTBp(p)STpBEAR	1539	2.6	c3
	VQHGAEp(p)STpEPR	1482	4.3	b10
	VAVpTPKSPK	1555	3.6	c1, c2
myristoylated a-lanine rich protein kinase C substrate	DKVTDGSGpSEK	1386	3.4	c2
	LSGpSPK	767	1.8	b7
	LSGpSPK	885.5	2.9	b11
Stathmin 1, leukemia associated phosphoprotein	VNGDpSPAAAEPAK	1337	3	b7
	DLpLESQK	1057	2.7	b7
	ApSGQAFELILSPR	1371	3.9	b7
MARCKS-like protein	SKp(p)SPDPLpSPRK	1510	4.2	b10, b11
	RApSGQAFELILSPR	1527	3.9	b11
SPR_MOUSE Scaffold attachment factor B2	LSLpSPRK	733	1.6	b7
	LSLpSPRK	896	2.9	b11
	AApTPESQEPQAK	1309	3	b7
Erythrocyte protein band 4.1, isoform 1	QDVTAEEAAGpSPAK	1355	4.3	b7
	GEVAPKETPK	1038	3.4	b10
	GEVAPKETPK	1165	2.3	c2
SPR_MOUSE Scaffold attachment factor B2	APTALpSPQDSK	1394	3.9	b7
Erythrocyte protein band 4.1, isoform 1	REAEQEVp(p)TPK	1571	3.9	c4

a) (p) refers to an uncertain phosphorylation site.
b) Xcorr lowest value 1.5, 2, 2.5 for charge states 1+, 2+, and 3+ respectively

Phosphorylated tyrosine does not lose phosphoric acid in the MS, i.e. it cannot be detected by neutral loss. Instead, the MS² spectra was studied. The sequence information is less in MS² compared to MS³ because of the lower collision energy. Despite this, some probable tyrosine-phosphorylated peptides were identified. A MS² spectra from a tyrosine-phosphorylated peptide is shown in Fig. 5.

The total number of proteins identified in the 31 fractions were approximately 5000, of which about 40 proteins (50 peptides) were found to be phosphorylated. Considering that about 10% of the proteins are phosphorylated, 500 phosphopeptides could theoretically be identified. For more identifications, prefractionation of the proteins should be performed before 2D LC on the peptides. Compared with the 1D LC method 10 times more phosphorylated peptides were found.

Conclusion

At least 50 phosphorylated peptides originating from 40 proteins were detected in the brain tissue with high confidence, many of the phosphorylation sites had not earlier been reported in the literature.

Ten times more phosphorylated peptides were found using 2D LC compared to 1D LC and with higher confidence.

Prefractionation of the proteins will be performed to increase the number of identifications.