

# Screening Medical Marijuana for Pesticides by GC-MS/MS

Di Wang, Laura Chambers, William Lipps, Zhuangzhi ‘Max’ Wang, Robert Clifford; Shimadzu Scientific Instruments, Columbia, MD, USA

## Introduction

Since medical marijuana (MM) was legalized in California in 1996, 23 states and Washington, D.C. have passed laws allowing its use for a variety of medical conditions. From a consumer safety point-of-view, quantitation of the pesticide residues in MM products has begun to attract wide interest. There are several problems associated with analysis of pesticide residues in MM. First and foremost, there are very few regulatory guidelines established to define which pesticides to include or what the detection limits should be, and secondly the matrix is very complex with significant interferences. Finally, sample load is growing exponentially, so the chosen method must be quick and easy to perform. Trace level pesticide analysis in complex food matrices have been done for many years with similar challenges, thus many of the analytical protocols emerging for the MM matrix are based on these well-established techniques.

Triple-quadrupole GC-MS/MS operated in MRM mode provides significant sensitivity and selectivity, but method development can be expensive and time consuming. This poster describes streamlined method development process for analysis of pesticide residues in MM using a QuEChERS sample preparation method, followed by GC-MS/MS detection and quantitation.

Note: Because medical marijuana has not been legalized in the state where the test lab is located, hops were used as the matrix in this application, as it is closely related to marijuana.

## Experimental

**Compound List**  
For this study 34 pesticides were selected for analysis based on the types of pesticides that are commonly used in MM production. The list includes several different compound classes (Table 1).

Organonitrogen Compounds	Synthetic Pyrethroid Compounds	Organophosphorus Compounds
Bupirimate	Bifenthrin	Chlorpyrifos
Etifenprox	Permethrin	Diazinon
Endosulfate (Terrazole)	Cyfluthrin	Malathion
Fenarimol	Deltamethrin	Mevinphos (Phosdrin)
Flutriafol	Flucythrinate	Phosalone
MKG-264	Lambda-cyhalothrin	Pirimiphos methyl
Myclobutanil	Tefluthrin	Carbamates and others
Paclobutrazol	Transfluthrin	2-Phenylphenol
Penconazole	Organochlorines compounds	Vinclozolin
Tebuconazole (Folicur)	Dichlorvos (DDVP)	
Terbutylazine	Endosulfan sulfate	
Triadimefon	gamma-BHC (Lindane)	
Triadimenol (Baytan)	p,p'-DDT	

Table 1: Selected Pesticide Compound Classes Included Organonitrogens, Synthetic Pyrethroids, Organochlorines, Organophosphates, and Carbamates

**Method Development**  
The most difficult part of any triple quadrupole method development process, is determination and optimization of the Multiple Reaction Monitoring (MRM) transitions and collision energies (CE). For this study, the Shimadzu Smart Pesticide Database was used as the foundation for creating the MRM analysis method. The Smart Pesticide Database includes up to six fully optimized MRM transitions and CEs for 479 pesticides and Retention Indices (RI) for accurately predicting compound retention times. The transitions and CEs in the database were optimized using the Shimadzu GCMS-TQ8040 triple quadrupole GC-MS/MS. Figure 1 shows a portion of the Smart Pesticide Database and the method, compound, and transition information.

Serial	Type	Acq. Mode	Method No.	Compound Name (S)	Ret. Index	CE	Transition	CE	Ratio
1	Target	MRM	1	Acetaminophen	115.1	100.0	Ref1	115.1	100.0
2	Target	MRM	1	Aspirin	180.1	100.0	Ref1	180.1	100.0
3	Target	MRM	1	Chlorzoxazone	281.1	100.0	Ref1	281.1	100.0
4	Target	MRM	1	Cyfluthrin	241.1	100.0	Ref1	241.1	100.0
5	Target	MRM	1	Deltamethrin	228.1	100.0	Ref1	228.1	100.0
6	Target	MRM	1	Endosulfan	321.1	100.0	Ref1	321.1	100.0
7	Target	MRM	1	Flucythrinate	228.1	100.0	Ref1	228.1	100.0
8	Target	MRM	1	Gamma-BHC	232.1	100.0	Ref1	232.1	100.0
9	Target	MRM	1	Malathion	102.1	100.0	Ref1	102.1	100.0
10	Target	MRM	1	Permethrin	259.1	100.0	Ref1	259.1	100.0
11	Target	MRM	1	Phosalone	228.1	100.0	Ref1	228.1	100.0
12	Target	MRM	1	Pyrimiphos methyl	228.1	100.0	Ref1	228.1	100.0

Figure 1: Example of Information Found in the Smart Pesticide Database Used to Create an MRM Analysis Method

A few of the target pesticides were not included in the Smart Pesticide Database. For these compounds, the MRM Optimization Tool was used to automatically determine the optimized MRM transitions and collision energies (CE). Once determined, the new transitions are added to the Smart Pesticide Database. Figure 2 shows the graphic results from the MRM Optimization Tool, with 6 transitions for two of the pesticides.




Figure 2: Optimized Transitions for Two Pesticides Using the MRM Optimization Tool

After adding the optimized transitions for the new pesticides to the existing Smart Pesticide Database, the MRM analysis method was created automatically. The program uses pesticide RIs in the database to accurately predict retention times for the target compounds. The Smart MRM function automatically adjusts Loop, Event, and Dwell times to optimize sensitivity for all compounds in the list simultaneously. Flexible MS events can create optimized methods with 400+ compounds. Used together, the Smart Pesticide Database and MRM Optimization Tool shortened the method development time from hours to just a few minutes.




Figure 3: The MRM Analysis Method is Created Automatically and Optimized for Sensitivity




Figure 4: Workflow using the MRM Optimization Tool


Information used to create the analysis method is shown in Table 2. It includes a compound table, retention indices and retention times, one transition with optimized CE for quantitation, and two reference transitions. Area ratios are also empirically determined, and can be used as part of the laboratory QA/QC program.

Serial	Compound Name	Ret. Index	Ret. Time	Type	ms	CE	Ratio	Type	ms	CE	Ratio	Type	ms	CE	Ratio			
1	Dichlorvos	1252	4.345	T	100.00	-79.00	7	100.00	Ref1	185.00	-93.10	13	44.15	Ref2	219.95	-185.00	15	10.19
2	Mevinphos	1427	5.642	T	127.05	-109.00	11	100.00	Ref1	192.05	-127.00	13	47.84	Ref2	127.05	-95.00	15	35.24
3	Endosulfate	1449	5.991	T	210.95	-183.00	11	100.00	Ref1	182.95	-140.00	15	96.56	Ref2	210.95	-180.00	23	91.67
4	2-Phenylphenol	1533	6.853	T	169.10	-141.10	13	100.00	Ref1	169.10	-112.10	25	91.99	Ref2	170.10	-141.10	23	96.59
5	Lindane	1779	8.660	T	180.95	-145.00	15	100.00	Ref1	218.90	-183.00	9	66.47	Ref2	218.90	-145.00	19	33.81
6	Terbutylazine	1782	8.694	T	229.10	-173.10	7	100.00	Ref1	214.10	-71.10	19	78.34	Ref2	214.10	-132.10	9	59.36
7	Diazinon	1790	8.766	T	304.10	-179.20	13	100.00	Ref1	240.05	-152.10	7	61.75	Ref2	240.05	-107.10	17	61.54
8	Tefluthrin	1816	9.002	T	177.05	-127.10	17	100.00	Ref1	177.05	-137.10	17	53.83	Ref2	197.05	-141.10	13	51.26
9	Vinclozolin	1894	9.730	T	212.00	-172.00	15	100.00	Ref1	212.00	-145.00	23	80.05	Ref2	285.00	-212.00	15	71.13
10	Transfluthrin	1903	9.815	T	163.05	-127.10	7	100.00	Ref1	163.05	-91.10	15	82.75	Ref2	163.05	-143.00	17	75.80
11	Metolaxyl	1915	9.926	T	234.10	-146.10	19	100.00	Ref1	234.10	-174.10	11	75.22	Ref2	249.15	-190.20	9	64.50
12	Permethrin methyl	1941	10.167	T	200.10	-125.10	23	100.00	Ref1	200.10	-233.10	11	51.89	Ref2	270.05	-125.00	17	54.23
13	Malathion	1964	10.377	T	127.05	-99.10	7	100.00	Ref1	173.10	-99.10	13	66.84	Ref2	173.10	-127.10	7	64.75
14	Chlorpyrifos	1980	10.529	T	313.95	-257.80	19	100.00	Ref1	315.95	-259.90	19	74.59	Ref2	285.95	-257.90	9	47.29
15	Triadimefon	2003	10.738	T	208.05	-111.10	23	100.00	Ref1	208.05	-127.10	15	89.54	Ref2	210.05	-143.10	9	43.85
16	MKG-264	2009	10.960	T	164.10	-93.10	13	100.00	Ref1	164.10	-96.10	13	65.56	Ref2	164.10	-90.10	23	55.12
17	Permethrin	2063	11.253	T	248.10	-157.10	25	100.00	Ref1	199.00	-123.10	19	50.14	Ref2	248.10	-192.10	15	45.77
18	Triadimenol	2092	11.541	T	168.15	-70.00	9	100.00	Ref1	128.00	-65.10	23	38.43	Ref2	112.05	-58.10	11	27.68
19	Paclobutrazol	2132	11.899	T	236.05	-125.10	11	100.00	Ref1	236.05	-167.10	9	37.10	Ref2	238.05	-127.10	11	32.47
20	Flutriafol	2135	12.104	T	235.00	-95.10	13	100.00	Ref1	219.05	-123.10	15	65.60	Ref2	123.05	-75.00	23	53.38
21	Myclobutanil	2209	12.952	T	179.05	-125.10	15	100.00	Ref1	179.05	-125.10	9	55.34	Ref2	179.05	-90.10	29	36.16
22	Buprimate	2204	12.535	T	273.10	-108.10	15	100.00	Ref1	273.10	-193.10	7	67.72	Ref2	193.15	-81.10	23	74.80
23	Endosulfan sulfate	2360	13.865	T	271.80	-236.80	21	100.00	Ref1	271.80	-234.90	17	22.20	Ref2	271.80	-141.00	13	22.31
24	p,p'-DDT	2367	13.919	T	225.00	-145.20	20	100.00	Ref1	237.00	-165.20	23	64.85	Ref2	235.00	-199.10	17	13.84
25	2-Chlorobenzonitrile	2399	14.184	T	250.10	-125.10	19	100.00	Ref1	250.10	-70.10	9	40.63	Ref2	252.10	-127.10	23	55.38
26	Bifenthrin	2471	14.767	T	181.15	-166.10	13	100.00	Ref1	181.15	-165.10	27	90.00	Ref2	166.10	-164.20	29	4.99
27	Phosalone	2556	15.432	T	182.05	-111.00	15	100.00	Ref1	182.05	-75.10	27	53.27	Ref2	182.05	-138.00	9	38.67
28	Lambda-Cyhalothrin	2597	15.748	T	197.05	-141.10	11	100.00	Ref1	208.10	-181.10	7	97.01	Ref2	197.05	-161.10	7	54.52
29	Permethrin	2621	16.067	T	251.00	-111.10	15	100.00	Ref1	251.00	-111.10	29	42.14	Ref2	300.05	-191.00	13	34.40
30	Permethrin	2706	16.562	T	183.00	-153.10	15	100.00	Ref1	183.00	-168.10	15	107.11	Ref2	163.00	-127.10	7	109.13
31	Cyfluthrin	2793	17.202	T	226.05	-206.10	13	100.00	Ref1	199.10	-170.10	25	70.95	Ref2	206.05	-151.10	19	64.85
32	Bifenthrin	2870	17.812	T	163.15	-135.10	11	100.00	Ref1	163.15	-107.10	17	89.29	Ref2	376.20	-163.20	11	5.79
33	Flucythrinate	2876	17.860	T	199.10	-157.10	9	100.00	Ref1	199.10	-107.10	23	94.17	Ref2	235.10	-119.10	19	18.37
34	Dichlorfenthrin	3061	19.609	T	252.90	-93.10	19	100.00	Ref1	181.10	-152.10	23	87.40	Ref2	252.90	-172.00	7	56.01

Gas Chromatograph	GC-2010 Plus
Injection Port	250 °C 1 µL splitless injection, 1 minute sampling time
Column	SH-Rxi-5Sil MS, 30 m x 0.25 mm x 0.25 µm film Helium carrier gas Constant Linear Velocity mode, 40.0 cm/second
Oven Program	85 °C (hold 1 minute) 25 °C/minute to 160 °C 10 °C/minute to 240 °C 10 °C/minute to 290 °C (hold 3 minutes)
Transfer Line	300 °C
Mass Spectrometer	GCMS-TQ8040
Acquisition Mode	MRM
Ion Source	230 °C Electron ionization mode, 70 eV
Collision Gas	Argon, 200 kPa
MRM Loop Time	Optimized with Smart MRM

Table 3: Optimized Instrument Conditions for Analysis of Pesticides in Hops (MM) Samples using the Shimadzu GCMS-TQ8040

### Sample Preparation - QuEChERS



### Calibration

A 5-point calibration curve was generated for all 34 target pesticides, covering the range from 1 to 200 parts-per-billion (ppb) (Figure 5). Figure 6 shows the overlaid MRM chromatograms from three transitions for two of the pesticides in the 1-ppb calibration standard.




Figure 5: Exponential Calibration Curves for Two Pesticides, 1 to 200 ppb




Figure 6: Example of Overlaid MRM Chromatograms For Two Pesticides in the 1-ppb Calibration Standard

### Sample Repeatability

Two different hops samples were processed using the QuEChERS procedure. The extracts were spiked with the pesticide mix at 25 ppb and analyzed in triplicate using the optimized MRM method. Chromatograms in Figure 7 illustrate how the MRM technique can be used to select the target compound from a complex matrix background, and produce reliable, reproducible results at low concentrations.




Figure 7: MRM Chromatograms of Two Hops Samples Spiked at 25 ppb and Analyzed in Triplicate

## Summary and Conclusion

The data presented illustrate how a triple quadrupole GC-MS/MS operated in the MRM mode, can be used to analyze for trace-level pesticide residues in complex plant matrices such as medical marijuana. The matrix was extracted using a QuEChERS kit, and interferences removed using an SPE cartridge. The resulting extracts were analyzed in triplicate using MRM transitions provided in the Smart Pesticide Database or individually optimized using the MRM Optimization Tool, with repeatability of 6% or better. The MRM method was fully optimized in just a few minutes, target compounds were selectively identified against the co-eluting matrix interferences, and quantitated at the parts-per-billion range.

## Acknowledgement

The authors wish to acknowledge Restek Corporation, Bellefonte, PA for useful discussions and advice regarding column selection, standards, and providing the QuEChERS materials used in the development of this method.