

MOHAWK Semi-automated tube picker

FEATURES

• Fast and reliable tube picking from 96 position racks

6

6

- Works with frozen and thawed tubes
- Create picklists in excel, csv, text, xml or json
- A 2D barcode scanner for tube confirmation and for generating picklists directly from the tube barcodes
- 1D scanner for rack barcode included as standard
- DataPaq[™]software with free lifetime upgrades

Easy to set up and Easy to Use

NA 14 14 14

The Mohawk works out of the box, and needs no set up or calibration. Instrument control is through user friendly software which provides an intuitive interface for easy operation.

Tubes can be picked from a single rack or a picklist can be set up to select tubes across multiple racks. Pick list operations can process racks from the linear barcode or manually select wells to pick.

SPECIFICATIONS

- Dimensions = 300mm x 250mm x 180mm (WxDxH)
- Total Weight = 5.7kg
- Power adaptor input = 100-240VAC
- Power adaptor output = 12VDC, 5A

- Picklist import file formats are excel, text, csv, json and xml
- Operating system on host PC = Windows 7, 8, 10
- Compatible with 96 position SBS rack

INCLUDED WITH THE MOHAWK

Mohawk tube picker instrument Power adaptor, power cable Mohawk software and license 1D rack code reader and bracket User guide

		HIERBER P
OPERATION		0
AUTOMATIC MODE		0
Designed to allow for operation of the Mohawk with	C	Zah
minimal interaction on control PC.		
Worklist File		
A file that describes the tubes that are to be selected,		
there are multiple formats available.		0

Ordering information

Code	Description	Quantity
ZTS-MHK	Ziath Mohawk Tube Picker	1

Ziath Ltd, Unit 2A Solopark Trading Estate, Pampisford, Cambridge, CB22 3HB, UK Ziath LLC, 8910 University Center Lane Aventine, Suite 400, San Diego, CA 92122 Tel UK: +44 (0) 1223 855021 Tel N.A.: (858) 880-6920 Web: www.ziath.com Email UK and N.A.: info@ziath.com

Ziath Japan Co., Ltd. BIZMARKS 603, Nihonbashi-Kayabacho, 8-2, Nihonbashi Koamicho, Chuo-ku, Tokyo 103-0016 Japan

Tel JPN: +81 3 6403 3359 Email JPN: info@ziath.co.jp