

Production of Prostanoids in the Plant Arabidopsis thaliana

Maged E. Mohamed and Colin M. Lazarus

School of Biological Sciences, University of Bristol, Woodland Road, Bristol BS8 1UG, United Kingdom, E-mail: m.e.mohammed@bristol.ac.uk



A. Introduction

Δ⁸ Pathway

Linoleic acid

(18:2 Δ^{9,12})

Eicosadienoic acid

(20-2 A11,14)

a son

Dihomo-y-linolenic

acid (20:3 Δ8,11,14)

Arachidonic acid

(20:4 ∆^{5,8,11,14})

Prostaglandin H₂

Prostaglandin F2

Prostaglandin L

TYAS

PGFS

A⁹ elongase

∆⁸ desaturase

⁵ desaturase

12 10 - 4

PGHS

Mouse PGHS cDNAs

amplification

Transformation

(0-

PCP

Yeast

PGHS-1

PGHS-2

• Prostanoids (prostaglandins (PGs) and thromboxanes) are lipid mediators derived from 20-carbon polyunsaturated fatty acids, which are found in nearly all mammalian tissues and fluids. They contribute to a number of physiological, pharmacological and pathological processes, such as inflammation, thrombosis, and gastro-intestinal tract secretions. • Transgenic plants provide cheap platforms for the production of high quality products in high yield. Their use eliminates the risk of contamination with human or animal pathogens such as HIV, hepatitis etc.

We previously achieved the production of di-homo-y-linolenic acid (DGLA), arachidonic acid (AA) and eicosapentaenoic acid (EPA) in Arabidopsis thaliana (Oi et al., 2004) via, so called, Δ^8 pathway by sequential transformation with a Λ^9 elongase from *Isochrysis galbana* (Oi *et al.*, 2002), a Λ^8 desaturase from *Euglena gracilis* (Wallis and Browse, 1999) and a Δ^5 desaturase from *Mortierella alpina* (Michaelson *et al.*, 1998). These 20-carbon fatty acids are the necessary feedstuffs for PG production via various prostanoid synthases.

B. Testing the activity of prostaglandin H synthase (PGHS) in yeast

• Mouse PGHS coding regions were first amplified from cDNA of the two PGHS isoforms, PGHS-1 and PGHS-2.


• The amplified isoforms were cloned in a yeast expression vector. Yeast cultures supplemented with DGLA and AA were induced and repressed by galactose and glucose respectively.

- The cells and cultures were assaved by enzyme immunoassay (EIA) to compare the activities of the different PGHS isoforms. • Results:
- PGHS genes were active and PGs were produced.
- PGHS-2 showed low activity relative to PGHS-1.

References

Qi B., Fraser T., Mugford S.,

- Most PG was found in the culture medium rather than in the cells.
- PGHS-1 showed the greatest activity with both substrates, while PHGS-2 displayed relatively good activity with DGLA as substrate (69 % of the PGHS-1 level) but very little activity above background with AA (11 % of the PGHS-1



pdfMachine

A pdf writer that produces quality PDF files with ease! Produce guality PDF files in seconds and preserve the integrity of your original documents. Compatible across Qi B., Beaudoinb F., Fraser T.
Wallis J. G.and Browse J. (199) Michaelson I V Lazanis C nearly all Windows platforms, simply open the document you want to convert, click "print", select the

"Broadgun pdfMachine printer" and that's it! Get yours now!

DGLA-producing


Δ thaliana

Agrobact

Selection

hygromyg

- C. Production of Prostaglandin Series 1 in A. thaliana
- DGLA-producing plants were transformed with PGHS-1 and PGHS-2.
- Transformed seeds were germinated on hygromycin-containing medium.
- Healthy seedlings were tested for the presence of the PGHS, Δ^8 desaturase and Δ^9 elongase genes by PCR.
- Plants showing positive results were transferred to soil and allowed to grow on and segregate. Individual homozygous lines were analyzed for PG production by EIA.
- Results: PGs were produced in PGHS-1 and to a lesser extent in PGHS-2 transformants. The most active PGHS-1 and PGHS-2 lines (D1-3 and H1-3. respectively) were propagated and 6 individuals of each line were analysed by EIA for PG production. The results (see below) show that the use of PGHS-1 is



better than PGHS-2 for this work. These lines are now being grown on a large scale for seed production and for extraction and purification of the prostaglandin produced so that it can be fully analysed.

D. Production of other types of prostanoids

- Mouse prostaglandin I synthase (PGIS) and thromboxane A synthase (TxAS) coding regions were amplified from cDNA.
- Trypanosoma brucei prostaglandin F synthase (PGFS) was amplified from genomic DNA.
- The three amplified genes were cloned into a yeast expression vector in sense and antisense orientations with respect to the GAL promoter, and then expressed in yeast together with PGHS-1 using AA as a substrate.
- Results: The three genes were active and prostaglandin I₂, thromboxane A, and prostaglandin $F_{2,a}$ were produced in various amount (see opposite)
- DGLA-producing plants are now being transformed with a Δ^5 -desaturase ene (to produce AA and EPA) together with PGHS-1 and each of PGIS,

S in order to produce the pharmaceutically active I_2 and F_{2a} or thromboxane A_2 , respectively.

